

Cold plates

AHP-1200CP

Environments

- Bench top, laboratory, industrial

Features

- Direct contact cooling as much as 48°C below room temperature
- Weighs only 19lbs (8.6kg)
- Compact bench top unit
- No compressor, fluorocarbons or filters
- Virtually maintenance-free operation
- Stainless steel exterior housing
- Integral temperature controller option
- Mounts in any orientation

Cooling capacity

- 260 watts @ 0°C triangle T

Options

- Adaptable for TC-3300 control
- Models with integral TC-3300 available
- Adjustable feet for FA-1200 option

Includes

- Integral power supply
- Cold plate accessory tapped holes
- Rubber feet
- Power input cord set
- Machined surface



Product no	Model no	Part number	Notes	Voltage Vac	Current Amps	Operation Ambient °C
12320/20	AHP-1200CP	1-3090-0-000	Cool only no control	120	4.0	-10/+70
12320/21	AHP-1200CP	1-3050-0-000	Cool only for remote control	120	4.0	-10/+70
12320/22	AHP-1200CP	1-30D0-0-000	Cool only TC-3300	120	4.0	-10/+70
12320/23	AHP-1200CPHC	1-30501-1-000	Heat/Cool for remote control	120	4.0	-10/+70
12320/24	AHP-1200CPHC	1-30D01-1-000	Heat/Cool TC-3300	120	4.0	-10/+70
12320/25	AHP-1202CP	1-3092-0-000	Cool only no control	240	2.5	-10/+70
12320/26	AHP-1202CP	1-3052-0-000	Cool only for remote control	240	2.5	-10/+70
12320/27	AHP-1202CPHC	1-3052-1-000	Heat/Cool for remote control	240	2.5	-10/+70