


# Containers and Vials

## Specimen Containers, Sterile


Polystyrene Sterile containers, aseptically produced, with caps in metal or polypropylene as indicated below. Supplied with a choice of no label, a plain white label, or a printed label which has a space for name, no., ward, specimen, date and time. Universal containers have a conical base, with skirt for ease of standing up.



Product No	Capacity	Cap Type	Label Type	Pack Qty
10005/01	30 (univ)	PP	none	400
10005/02	30 (univ)	PP	plain	400
10005/03	30 (univ)	PP	printed	400
10005/04	60	PP	none	300
10005/05	60	PP	plain	300
10005/06	60	PP	printed	300
10005/07	60	metal	none	300
10005/08	60	metal	plain	300
10005/09	60	metal	printed	300
10005/10	100	metal	none	120
10005/11	100	metal	plain	120
10005/12	100	metal	printed	120
10005/13	150	metal	none	90
10005/14	150	metal	plain	90
10005/15	150	metal	printed	90
10005/16	250	metal	none	48
10005/17	250	metal	plain	48
10005/18	250	metal	printed	48

## Clarified Polypropylene


Greater physical durability & chemical resistance, combined with good optical quality



Product No	Capacity	Cap Type	Label Type	Pack Qty
11036/01	30 (univ)	PP	none	400
11036/02	30 (univ)	PP	plain	400
11036/03	30 (univ)	PP	printed	400
11036/04	60	PP	none	300
11036/05	60	PP	plain	300
11036/06	60	PP	printed	300
11036/07	60	metal	none	300
11036/08	60	metal	plain	300
11036/09	60	metal	printed	300
11036/10	100	metal	none	200
11036/11	100	metal	plain	200
11036/12	100	metal	printed	200
11036/13	150	metal	none	100
11036/14	150	metal	plain	100
11036/15	150	metal	printed	100
11036/16	250	metal	none	50
11036/17	250	metal	plain	50
11036/18	250	metal	printed	50

## Hinged-Top Containers

A range of leak-proof aseptically produced containers manufactured from translucent polypropylene. The captive cap allows one-handed operation, and the patented hinge mechanism holds the cap open during sampling or analysis. Containers are available with a pierceable diaphragm, allowing sampling without opening the cap.



Product No	Capacity	Colour	Pack Qty
11248/01	60	Transparent	650
11248/02	60	Blue	650
11248/04	100	Transparent	350
11248/05	100	Blue	350
11248/06	100	Transparent	350
11248/07	300	Transparent	240

