


# Containers and Vials

## Trident Vials

Made from neutral borosilicate glass. Supplied WITH caps, bagged separately.



### Rolled Rim, Squat Form

Product No	Approx Capacity-ml	Height-mm	Diameter-mm	Approx Pack Qty
10095/01	7	34	23	190
10095/02	10.5	46	23.25	190
10095/03	14	48	25.5	162
10095/04	21.25	66	25.5	162
10095/05	28.25	72	27	140

### Tall Form, Screw Neck

Product No	Approx Capacity-ml	Height-mm	Diameter-mm	Approx Pack Qty
10093/01	1.75	36	11	882
10093/02	3.5	46	12.5	666
10093/03	7.0	58	17	399
10093/04	10.5	67	17	399
10093/05	14.0	72	20	264
10093/06	21.25	86	21	231
10093/07	28.25	94	23	190

### Pathology Vial, Screw Neck

Product No	Approx Capacity-ml	Height-mm	Diameter-mm	Approx Pack Qty
10094/01	7	42	20	264
10094/02	14	58	23.25	190
10094/03	28.25	72	27.5	140