

## EC counting scale - A fast, accurate and user friendly counting machine

### Applications

- Weighing (3 units)
- Parts counting with average piece weight and auto optimisation
- Checkweighing
- 80 hour internal rechargeable battery with auto shut off

### Features

- Removable stainless steel weighing pan and durable plastic housing
- 3-window, backlit LCD for fast display of weight
- 9 function and 11 numeric keys plus dedicated TARE and SAMPLE keys
- Accumulator function for weights and parts count
- User programmable parts checking mode
- 10 APW (Average piece weight) memory locations
- Adjustable filtering levels and auto zero tracking
- Rubber levelling feet and front level bubble
- Supplied with protective cover
- Refer to later pages for the correct calibration weights

Product no	Capacity	Readability	Platform
<b>12206/50</b>	3kg	0.1g	300 x 230mm
<b>12206/51</b>	6kg	0.2g	300 x 230mm
<b>12206/52</b>	15kg	0.5g	300 x 230mm
<b>12206/53</b>	30kg	1g	300 x 230mm

### Accessories

Product no	Description
<b>12206/55</b>	Replacement mains adapter



## ESTM series

The multipurpose low profile bench scale

### Applications

- Weighing (4 units)
- Mains or battery operation

### Features

- Two button operation
- Stability indicator
- Auto-zero tracking
- Auto-off timer
- 3 way mounting bracket (column, wall and platform)
- Refer to later pages for correct calibration weights

Product no	Capacity	Readability	Platform
<b>12206/60</b>	6kg	2g	310 x 270mm
<b>12206/61</b>	30kg	10g	310 x 270mm
<b>12206/62</b>	50kg	20g	310 x 270mm
<b>12206/63</b>	50kg	20g	520 x 400mm
<b>12206/64</b>	100kg	50g	520 x 400mm
<b>12206/65</b>	200kg	100g	520 x 400mm

### Accessories

Product no	Description
<b>12206/67</b>	Replacement AC mains adapter

