

## Liquid Nitrogen Dewars

### Cryostor Series

- Stainless steel cylinders
- Castor mounted for easy manoeuvrability
- Dual decant facility on manifold
- Faster liquid and gas delivery rates
- Standard pressure building circuit with regulator, pressure gauge, relief valve and bursting disc
- Suitable for Liquid Nitrogen and Argon
- Analogue contents gauge indicator



Product No	Model	Capacity-Litres	Overall Height-mm	Diameter-mm	Weight (empty)-kg	Daily evaporation rate
11112/01	Crystor 30	30	800	461	61	3.2%
11112/02	Crystor 60	60	925	537	113	1.8%
11112/03	Crystor 90	90	1070	537	138	1.4%
11112/04	Crystor 120	120	1107	640	192	1.0%
11112/05	Crystor 180	180	1386	640	261	0.8%
11112/06	Crystor 240	240	1651	640	329	0.7%

### Cryostor Accessories

Product No	Description
11112/11	4 foot flexible hose
11112/12	5 foot flexible hose
11112/13	2 metre flexible hose
11112/14	1/2" BSP Phase separator

## Open Dewars

### OD Series

- Stainless steel vacuum insulated dewars
- Complete with lid and carrying handle, except (OD 0.5)
- Resistant to vibration, impact or breakage
- Easy to use
- Lightweight models 0.5 - 6 litre
- Heavy duty models 10-80 litre
- Suitable for use with cold traps, oil baths, water baths and in cryosurgery

Product No	Model	Capacity-Litres	Duty
11113/01	OD 0.5	0.5	Lightweight
11113/02	OD 1	1	Lightweight
11113/03	OD 2	2	Lightweight
11113/04	OD 3	3	Lightweight
11113/05	OD 5	5	Lightweight
11113/06	OD 6	6	Lightweight
11114/01	OD10	10	Heavy duty
11114/02	OD15	15	Heavy duty
11114/03	OD25	25	Heavy duty
11114/04	OD50	50	Heavy duty
11114/05	OD80	80	Heavy duty

