


## Cryosleeves



PVC clear plastic sleeves enclose cryocanes (Cat No 10475) for additional security during storage and transportation. Can be used repeatedly and are much more durable than conventional cardboard

Product No	Pack Qty
10476/01	100

## Racking Systems



A range of stainless steel racks for use with the 11120 range polypropylene boxes, suitable for use in chest or upright freezers if required. Models 11122/01 to 11122/06 have a front retaining bar to hold boxes in place, such as whilst being withdrawn from a chest freezer.

Product No	No of Shelves	H x W x D (mm)
11122/01	6	404 x 147 x 154
11122/02	7	470 x 147 x 154
11122/03	8	527 x 147 x 154
11122/04	8	555 x 147 x 154
11122/05	11	732 x 147 x 154
11122/06	12	785 x 147 x 154
11122/11	8	280 x 300 x 150

## Racking Systems, with Drawers

Product No	No of Drawers	Box Capacity	H x W x D-mm
11122/21	6	18	400x150x460
11122/22	6	24	400 x 150x612

## Storage Boxes, Ultra Low

A range of stackable polycarbonate boxes that may be used at temperatures down to -196°C and can withstand autoclaving at 121°C. Optional inserts will accept most 1ml and 2ml cryogenic vials and low volume microtubes making them particularly versatile. The lids of the 100 tube boxes have an embossed grid with numbers (1-100). All boxes are supplied in packs of five. Overall dimensions: 50 x 130 x 130mm. Supplied in packs of 5.

Product No	Lid	Insert	Tube Dia-mm	Spaces
11123/01	Gridded	White	12 - 12.5	100
11123/02	Gridded	Yellow	10 - 10.5	100
11123/03	Clear	Removable	12 - 12.5	81
11123/04	Clear	None	N/a	N/a

## Storage Boxes & Racks

A range of polypropylene boxes and racks that can be used independently or together to suit various needs. They withstand autoclaving at 121°C for 15 minutes and storage at -90°C. The box lids and racks include an alpha-numeric grid reference for ease of tube identification. These stackable boxes have 100 positions suitable for use with 0.5ml and 1.5/2.0ml microtubes as well as most manufactured brands of cryovials. The maximum tube height is 53.5mm. The boxes have highly translucent lids in a choice of colours, with natural bases, or all black (opaque) where total light exclusion is required. (The bases and lid sides are mat to accept marking). Each box measures 140 x 140 x 60mm. The 64-position racks are suitable for 0.5ml or 1.5ml microtubes (or other items of similar diameter). They can be used alone, or in the boxes to simplify tube handling in batches. When stood on the bench the tubes will generally be raised slightly to aid individual selection.

Product No	Colour
11120/01	Natural
11120/02	Red
11120/03	Yellow
11120/04	Blue
11120/05	Green
11120/06	All black

## Insert / Racks

### for 0.5ml Tubes

Product No	Colour
11121/01	Natural
11121/02	Red
11121/03	Yellow
11121/04	Blue
11121/05	Green
11121/06	Black

### for 1.5ml Tubes

Product No	Colour
11121/11	Natural
11121/12	Red
11121/13	Yellow
11121/14	Blue
11121/15	Green
11121/16	Black