

Scalpel Blade Remover

A single-use, safe and simple method for removing and disposing scalpel blades. Holding the handle, introduce the blade into the remover. Close it securely and remove the handle. Discard the remover with blade inside. Suitable for use with handles No. 3 and 4, and all blades except No.12. Sterile removers are individually packed.

Product no	Sterile	Pack Size
11129/01	No	100
11129/02	No	500
11129/03	yes	50

Scalpel Handles

Non-sterile, nickel alloy handles. Handle No's 3, 3L and 7 fit blades 6 to 16. Handle No4 and 4L fit blades 18 to 27.

Product no	Handle No
10201/01	3
10201/04	3L
10201/02	4
10201/05	4L
10201/03	7*

*Handle No. 7 formerly No. 5

Scalpel, Sterile Disposable

Polystyrene handled disposable scalpels in individual sleeves. Gamma, irradiated. Supplied in packs of 10. Other blade shapes available to order.

Product no	Blade Shape
10200/01	10
10200/02	10A
10200/03	11
10200/04	12
10200/05	15
10200/06	15A
10200/07	20
10200/09	22
10200/10	22A
10200/12	24

Post Mortem Knife

A die cast handle with screw-type fitting and VPI protected heavy duty carbon steel individually wrapped blades. (Supplied separately).



Product no	Description	Pack size
10204/01	Handle	1
10204/02	Blade (PM40)	10

Scalpel

Product no	Description
10202/01	A handle of impact polystyrene with one No.1 blade and one No. 2 blade.



Spare Blades (in packs of 50)

Product no	Blade shape
10203/01	1
10203/02	2
10203/03	3




Scissors Dissecting

Straight

Stainless steel scissors with straight points and open shanks. Supplied by the pair in the lengths indicated below.

Product no	Length (mm)
10334/01	112
10334/02	125
10334/03	150



Curved

Stainless steel scissors with curved points and open shanks. Supplied by the pair in the lengths indicated below.

Product no	Length (mm)
10335/01	112
10335/02	125
10335/03	150

