

Drying

Neildorn Glassware Dryer



Thermostatic control of heating elements.

Takes up to 60 items.

Can be either wall mounting or bench standing.

Fitted with two heat switch plus cold, a detachable mains lead 1.5 metres long complete with 13 Amp plug and 30 minute timer. 220/240 volt 50Hz

Product no 10436/01



Draining Rack, Plastic

A one-piece high impact polystyrene backplate (450 x 630mm) incorporates a drain channel which may be connected to a drain tube (provided). The 72 pegs measure 95mm long x 15mm diameter and are removable for cleaning. A wall-mounting kit is included, and an optional pack of eleven small diameter pegs (95mm x 6mm) is available for test tubes, etc.



Product no	Description
10525/01	Rack
10525/02	6mm pegs (pack of 11)

Drying/Warming cabinets

Designed for drying glassware or instruments, these cabinets offer a temperature range up to 50°C above ambient. (The DC/540 and DC/1000 cabinets are also available in 65°C above ambient). The cases are stove powder painted for durability and ease of cleaning, and the doors are of toughened glass. The shelves are of chromium plated grid wire and are adjustable on a ladder racking system.



Product no	Capacity (litres)	Internal H x W x D (mm)	External H x W x D (mm)	No of Shelves	Wattage	Convection	Variable Rate extraction	Door style
11152/01	113	460 x 728 x 341	559 x 730 x 360	3	750	Natural	No	2 x sliding
11152/02	226	567 x 926 x 430	666 x 928 x 449	3	1000	Natural	No	2 x sliding
11152/03	540	1550 x 581 x 600	1664 x 583 x 610	3	1750	Natural	No	1 x hinged
11152/04	540	1550 x 581 x 600	1814 x 583 x 610	3	1750	Fan assisted	No	1 x hinged
11152/05	540	1550 x 581 x 600	1924 x 583 x 610	3	1750	Natural	Yes	1 x hinged
11152/06	540	1550 x 581 x 600	1924 x 583 x 610	3	1750	Fan assisted	Yes	1 x hinged
11152/07	1000	1550 x 1108 x 575	1664 x 1110 x 585	4	2750	Natural	No	2 x hinged
11152/08	1000	1550 x 1108 x 575	1814 x 1110 x 585	4	2750	Fan assisted	No	2 x hinged
11152/09	1000	1550 x 1108 x 575	1924 x 1110 x 660	4	2750	Natural	Yes	2 x hinged
11152/10	1000	1550 x 1108 x 575	1924 x 1110 x 660	4	2750	Fan assisted	Yes	2 x hinged