

Extraction

TURBOSOG

All digestion systems have a water heat pump WSP which removes fumes quite efficiently within a fume cabinet. For removal of aggressive acid fumes outside a cabinet we recommend the use of the Turbosog-scrubber unit.

Turbosog-Scrubber Unit

The Turbosog centrifugal scrubber condenses and neutralises aggressive acid fumes. The preseparator works in two steps, separating and washing out these acid fumes. The removal of the fumes is very efficient with extremely low running costs. Little service is required as no activated carbon filters are used in the system.

Product no	Type	Description
12550/60	TUR/K	Turbosog - Scrubber



KJELTABS

Catalyst tablets.

Tin with 1000 tablets

Product no	Type	Description
12550/61	ST	3.5g K ₂ SO ₄ + 0.0035g Se
12550/62	S	5.0g K ₂ SO ₄ + 0.005g Se
12550/63	CX	5.0g K ₂ SO ₄ + 0.5g CuSO ₄ x 5H ₂ O
12550/64	CT	5.0g K ₂ SO ₄ , 0.15g CuSO ₄ x 5H ₂ O + 0.15G TiO ₂
12550/65	SQ	1.5g K ₂ SO ₄ + 0.0015g Se
12550/66	IB/61	5.0g K ₂ SO ₄ , CuSO ₄ x 5H ₂ O + Se (100:6:1 parts) acc. to Weininger
12550/67	CK	1.5g K ₂ SO ₄ + 0.4g CuSO ₄ x 5H ₂ O
12550/68	TCT	3.5g K ₂ SO ₄ + 0.105g CuSO ₄ x 5H ₂ O + 0.015G TiO ₂ x 5H ₂ O
12550/69	CQ	1.5g K ₂ SO ₄ + 0.15g CuSO ₄ x 5H ₂ O
12550/70		K ₂ SO ₄ + CuSO ₄ x 5H ₂ O
12550/71	KS	100 x K ₂ SO ₄ 1 x Se

Antifoam tablets

Tin with 1000 tablets

Product no	Type	Description
12550/72	AS	Antifoam tablets

