

Technical Data

Type	TT 625/TT 625M	TT 125/TT 125M	TT 440/TT 440M	TT 480/TT 480M
Product no	12550/20 & 21	12550/22 & 23	12550/24 & 25	12550/26 & 27
Nominal Voltage	230 V AC	230 V AC	230 V AC	230 V AC
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Nominal Wattage	1500W	1500W	1500W	1500W
Weight	21kg	22.5kg	22.5kg	21.5kg
Dimensions (w x d x h)	525 x 450 x 740mm	525 x 450 x 740mm	525 x 450 x 740mm	525 x 450 x 740mm
Temperature	750°C	750°C	750°C	750°C
Heating places	6	12	4	4
Size of tubes	250ml	250ml	400ml	800ml
Type	TT 100/TT100M	KB 8 S	KB 8 S-BS	KB 20 S
Product no	12550/28 & 29	12550/50	12550/51	12550/52
Nominal Voltage	230 V AC	230 V AC	230 V AC	230 V AC
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Nominal Wattage	1500W	1000W	1000W	2200W
Weight	22.5kg	16kg	16kg	26kg
Dimensions (w x d x h)	525 x 450 x 740mm	380 x 380 x 650mm	380 x 380 x 650mm	380 x 380 x 650mm
Temperature	750°C	430°C	430°C	430°C
Heating places	12	8	8	20
Size of tubes	100ml	250ml	400ml	250ml
Type	KB 40 S	KBL 8 S	KBL 8 S-BS	KBL 20 S
Product no	12550/53	12550/54	12550/55	12550/56
Nominal Voltage	230 V AC	230 V AC	230 V AC	230 V AC
Frequency	50/60Hz	50Hz	50Hz	50Hz
Nominal Wattage	2200W	1160W	1160W	2360W
Weight	26kg	29kg	29kg	39kg
Dimensions (w x d x h)	380 x 500 x 650mm	460 x 500 x 740mm	460 x 500 x 740mm	460 x 500 x 740mm
Temperature	430°C	430°C	430°C	430°C
Heating places	40	8	8	20
Size of tubes	100ml	250ml	400ml	250ml
Type	KBL 40 S	CSB 8 M	CSB 20 M	CSB 8 A
Product no	12550/57	12550/30	12550/31	12550/34
Nominal Voltage	230 V AC	230 V AC	230 V AC	230 V AC
Frequency	50Hz	50/60Hz	50/60Hz	50Hz
Nominal Wattage	2360W	1000W	2200W	1160W
Weight	39kg	15kg	24kg	31kg
Dimensions (w x d x h)	460 x 500 x 740mm	390 x 390 x 1150mm	390 x 520 x 1150mm	470 x 390 x 800mm
Temperature	430°C	430°C	430°C	430°C
Heating places	40	8	20	8
Size of tubes	100ml	250ml	250ml	250ml
Type	CSB 20 A	SMA 8 A	SMA 20 A	
Product no	12550/35	12550/36	12550/37	
Nominal Voltage	230 V AC	230 V AC	230 V AC	
Frequency	50Hz	50Hz	50Hz	
Nominal Wattage	2360W	1160W	2360W	
Weight	45kg	32kg	46kg	
Dimensions (w x d x h)	470 x 520 x 800mm	470 x 390 x 1000mm	470 x 520 x 1000mm	
Temperature	430°C	430°C	430°C	
Heating places	40	8	20	
Size of tubes	100ml	250ml	250ml	
Type	TR	TUR/K		
Product no	12550/58	12550/60		
Nominal Voltage	230 V AC	230 V AC		
Frequency	50/60Hz	50/60Hz		
Nominal Wattage	-	205W		
Weight	1.4kg	19kg		
Dimensions (w x d x h)	80 x 150 x 155mm	330 x 450 x 420mm		
Temperature	0 - 430°C	-		
Accuracy	<0.5%	-		
Programs	-	-		