


# Extraction Accessories

## Soxhlet Extraction Heaters


Electrothermal EME series soxhlet heaters are available in three or six positions for four flask sizes as indicated below. All units have a "power on" indicator light, and each position has a "heater on" indicator light and controller for differing temperature requirements up to 450°C. Supporting rods are supplied with all units, and each position has an earth screen. These units are powder coated for chemical resistance.



Product no	Flask Size (ml)	No of Recesses	Wattage
10565/01	100	3	180
10565/02	250	3	450
10565/03	500	3	600
10565/04	1000	3	900
10565/05	100	6	360
10565/06	250	6	900
10565/07	500	6	1200
10565/08	1000	6	1800

## Soxhlet Extraction Heaters, Stirred

Similar to Product No 10565 above but with integral stirring to speed extractions. Stirring speed is continuously variable from 50-1000rpm via a single controller on the 3-place heaters or two controllers on 6-place heaters, each controlling three positions. Supplied for the flask capacities indicated below.



Product no	Flask Size (ml)	No of Recesses	Wattage
10566/01	100	3	280
10566/02	250	3	550
10566/03	500	3	700
10566/04	1000	3	1000
10566/05	100	6	460
10566/06	250	6	1000
10566/07	500	6	1300
10566/08	1000	6	1900

## Extractors, Soxhlet, BS2071

Manufactured from borosilicate glass.



Product no	Nominal Capacity (Syphoning Volume) (ml)	Socket Size	Cone Size	Approx. Overall Height (mm)	Recom'd Thimble Size (mm)
10670/01	6	19/26	14/23	150	10 x 50
10670/02	20	24/29	19/26	175	18 x 65
10670/03	40	29/32	24/29	200	22 x 80
10670/04	60	34/35	24/29	205	28 x 80
10670/05	60	34/35	34/35	205	28 x 80
10670/06	100	34/35	34/35	240	28 x 120
10670/07	100	40/38	24/29	220	33 x 100
10670/08	200	50/42	24/29	250	41 x 123

\*NB Height is measured from wide end of cone.

## Extractors, Soxhlet, BS2071

Similar to Product no 10670 above but with a flat flange joint at the top of the body for easy loading. Suitable lids and flange clips should be ordered separately, as required.

Product no	Nominal Capacity (Syphoning Volume)	Flange Size	Cone Size
10671/01	600 ml	75 mm	34/35
10671/02	2000 ml	100 mm	34/35

## Extraction Accessories

Borosilicate Glass Lids:



Product no	Flange Bore (mm)	Centre Socket	Parallel Socket	10° Side Sockets
10682/02	75	34/35	-	-
10682/10	100	34/35	-	-
10682/11	100	34/35	19/26	2 x 19/26

Flange Clips (Wire):

Product no	Flange Bore (mm)
10683/01	75
10683/02	100

