

Extraction Accessories

These plain extractor bodies can be used with Soxhlet liners or extraction thimbles. Conical glass spirals to support the sample or thimble should be ordered separately, if required.



Product no	Nominal Capacity (ml)	Socket Size (mm)	Cone Size	Approx Overall Height (mm)
10672/01	60	34/35	24/29	205
10672/02	100	45/40	24/29	225
10672/03	200	55/44	24/29	250

Glass Spirals:

Product no	To Suit Extractor (ml)
11161/01	60
11161/02	100 & 200

Soxhlet Liners

Borosilicate glass syphon assembly which, when used with plain body extractors (Cat No 10672 above) form a jacketed Soxhlet.

Product no	Nominal Capacity (Siphoning Volume)	To fit Body
10673/01	20ml	10672/01
10673/02	50ml	10672/03
10673/03	100ml	10672/03

Kjeldahl Catalyst Tablets

A wide range of tableted reagents is available to order including some 30 Kjeldahl tablets for specific applications. Most are made to order (please enquire), however, below are listed the most commonly used. All are supplied in packs of 1000.

Product No	Description	Tablet contents
10740/01	Sodium sulphate	1.0g Sodium sulphate
10740/02	Selenium	1.0g Sodium sulphate & 0.05g Selenium
10740/03	Copper	1.0g Sodium sulphate & 0.1g Copper II sulphate 5-hydrate
10740/04	Mercury	1.0g Sodium sulphate & 0.108g Mercuric oxide
10740/05	Titanium	1.0g Sodium sulphate & 0.03g Copper II sulphate 5-hydrate & 0.03g Titanium dioxide

Extraction Accessories

Thimbles of the highest quality, of pure cellulose. They are seamless for maximum strength.

These are suitable for general use such as for the extraction of organic compounds from reaction mixtures, foodstuffs, natural materials, analysis of lacquer and binder in paints, dust and tar determinations in tobacco, etc. Having a thickness of 2mm, they are similar to Whatman double thickness thimbles.

Product no	Inner Dia x Length (mm)	Product no	Inner Dia x Length
12590/01	9x40	12590/30	33x90
12590/02	10x50	12590/31	33x94
12590/03	16x50	12590/32	33x100
12590/04	16x100	12590/33	33x118
12590/05	19x90	12590/34	33x120
12590/06	22x80	12590/35	33x130
12590/07	25x60	12590/36	33x205
12590/08	25x70	12590/37	34x130
12590/09	25x80	12590/38	35x120
12590/10	25x90	12590/39	35x150
12590/11	25x100	12590/40	40x85
12590/12	26x60	12590/41	40x150
12590/13	27x80	12590/42	41x123
12590/14	27x100	12590/43	41x130
12590/15	28x60	12590/44	43x123
12590/16	28x80	12590/45	48x145
12590/17	28x100	12590/46	48x200
12590/18	28x120	12590/47	44x230
12590/19	30x80	12590/48	53x315
12590/20	30x90	12590/49	60x275
12590/21	30x100	12590/50	75x250
12590/22	33x60	12590/51	75x330
12590/23	33x80		

Economy

Similar to above, but only half the thickness (approximately 1.0mm).

Product no	Inner Dia x Length (mm)
12590/60	19x90
12590/61	26x60
12590/70	28x80
12590/71	30x80