

Filtration

Filter Papers - Grade No 50

Typical Applications Qualitative analysis of acid and alkaline solutions, fibre retention papers, vacuum & pressure filtration

Properties Hardened, slow, wet strength

Product No	Diameter-mm	Pack Qty
12602/90	55	100
12602/91	90	100
12602/92	110	100
12602/93	125	100
12602/94	150	100
12602/95	185	100
12602/96	240	100

Filter Papers - Grade No 52

Typical Applications Qualitative analysis of acid and alkaline solutions, fibre retention papers, vacuum & pressure filtration

Properties Hardened, medium fast, wet strength

Product No	Diameter-mm	Pack Qty
12602/70	55	100
12602/71	110	100
12602/72	150	100
12602/73	185	100
12602/74	580 x 580	500
12602/80	110	100
12602/81	150	100
12602/82	185	100
12602/83	240	100
12602/84	270	100
12602/85	320	100
12602/86	385	100

Filter Papers - Grade No 91

Typical Applications Filtering fine contaminations; use as protective paper for filter press cloths, cocoa butter, edible oils

Properties Medium to slow, wet strength

Product No	Diameter-mm	Pack Qty
12602/30	110	100
12602/31	150	100
12602/32	185	100
12602/33	580 x 580	500

Filter Papers - Grade No 54

Typical Applications Qualitative analysis of acid and alkaline solutions, fibre retention papers, vacuum & pressure filtration

Properties Hardened, fast, wet strength

Product No	Diameter-mm	Pack Qty
12602/50	55	100
12602/51	90	100
12602/52	110	100
12602/53	125	100
12602/54	150	100
12602/55	185	100
12602/56	25.5	210m roll
12602/60	110	100

12602/61	125	100
12602/62	150	100
12602/63	185	100
12602/64	240	100

12602/65	270	100
12602/66	32	100
12602/67	385	100
12602/68	500	100

Filter Papers - Grade No 113

Typical Applications Agar-agar, alginate, essential oils, colloids, gelatine, resin solutions, lacquers, syrupy juices, syrups, sweetened spirits.

Properties Very fast, high wet strength, thin

Product No	Diameter-mm	Pack Qty
12600/01	110	100
12600/02	150	100
12600/03	185	100
12600/04	240	100
12600/05	270	100

Filter Papers - Grade No 541

Typical Applications Quantitative filtrations of acid and alkaline solutions; pressure and vacuum filtrations

Properties Hardened, fast, wet strength

Product No	Diameter-mm	Pack Qty
12602/40	70	100
12602/41	90	100
12602/42	110	100
12602/43	125	100
12602/44	150	100
12602/45	185	100
12602/46	240	100