

Postlip Filter Papers

Smooth cellulose Thin White filter papers for general purpose and education use. The following have a weight of 70g/m² and a retention rating of 6µm. Other grades are available.

Product No	Diameter/ Size (mm)	Pack Quantity
11162/01	95	100
11162/02	115	100
11162/03	125	100
11162/04	150	100
11162/05	190	100
11162/06	230	100
11162/11	500 x 500	500

Seed Test Paper

Product No	Grade	Size (mm)	Description	Pack Quantity
12603/01	0858	580 x 110	Wrapping strips for pleated strips	500
12603/02	3014	2000 x 110	Pleated white strips for germination testing of beet seed after Neeb	1000
12603/03	3236	2000 x 110	Pleated grey strips for germination testing of seeds with fine germination tubes after Neeb	1000



Glass Fibre Filter Paper

All Schleicher & Schuell glassfibre filters are made of 100% borosilicate glass for high chemical resistance to a wide range of organic and inorganic solutions. The fibres are of mixed lengths and diameters for a depth filter structure and a large internal filter area. This gives these filters an outstanding pollutant-capturing capacity and low filter resistance.

Grade No GF/A

Product No	Diameter/ Size (mm)	Pack Quantity
12603/10	25	200
12603/11	47	200
12603/12	50	200
12603/13	55	100
12603/14	70	100
12603/15	90	100
12603/16	110	100
12603/17	125	100
12603/18	150	100
12603/19	625 x 660	100

Paper - Grade No GF/C

Typical application

Analysis of particulate or sludge fractions in waste and potable water; radio-immunoassays

Product No	Diameter/ Size (mm)	Pack Quantity
12603/20	25	200
12603/21	47	200
12603/22	50	200
12603/23	55	100
12603/24	70	100
12603/25	90	100
12603/26	110	100
12603/27	125	100
12603/28	625 x 660	100

Chromatography Papers

Grade No 3MM Chr

Product No	Size (mm)	Pack
12603/30	460 x 570	100
12603/31	580 x 600	100