

Balances

BW Series - Washdown IP65 Compact Bench Scale

The Ohaus BW series Compact Bench Scale is a cost-effective portioning, packing and general weighing scale designed to provide fast and accurate weighing performance in dry and washdown environments. With an internal sealed rechargeable lead-acid battery, IP65 housing and high-speed display update, the Ohaus BW Series offers durability and dependability in an affordable design.

Applications

- Weighing (2 units)
- Over-Under checkweighing function
- Percentage weighing
- Mains adaptor or internal rechargeable lead-acid battery

Features

- 1 red LED weight display: 6-digit, 7-segment, 14mm high characters, adjustable brightness level
- Fast stabilisation ≤ 2 seconds
- 4 raised tactile mechanical buttons: On/Off, Check, Tare, Zero
- Durable high-profile ABS housing with stainless steel pan (pan is dishwasher safe)
- 4 adjustable rubber levelling feet with integral front-panel levelling indicator
- 50-hour internal rechargeable battery with power-saving Auto-shut off feature
- Battery level indication
- IP65 dust and water protection
- Max. displayed resolution from 1:6000-7500 non approved, 1-3000 approved



Model no	Product no	Capacity	Readability (non approved)	Readability (approved)	Platform (mm)
BW1.5	12210/01	1.5 kg	0.2g	0.5g	230 x 180
BW1.5M*	12210/02	1.5 kg	0.2g	0.5g	230 x 180
BW3	12210/03	3 kg	0.5g	1g	230 x 180
BW3M*	12210/04	3 kg	0.5g	1g	230 x 180
BW6	12210/05	6 kg	1g	2g	230 x 180
BW6M*	12210/06	6 kg	1g	2g	230 x 180
BW15	12210/07	15 kg	2g	5g	230 x 180
BW15M*	12210/08	15 kg	2g	5g	230 x 180

Accessories

Product no	Description
12210/09	Large platform

*EC verified

*1 Postcode of final destination must be supplied at time of order