

Bottle-Top Filter Units, Sterile

Similar to the above units, but WITHOUT receivers. Suitable media bottles will be found in the Duran range and must match the cap size. Details of the membrane properties will be found previously. Supplied in packs of 12.



SFCA Membrane (Yellow Collar)

Product No	Reservoir Vol. (ml)	Pore size (µm)	Membrane Dia. (mm)	Cap size
10550/01	150	0.20	50	GL32
10550/02	150	0.45	50	GL32
10550/03	150	0.20	50	GL45
10550/04	50	0.45	50	GL45
10550/05	500	0.20	75	GL32
10550/06	500	0.45	75	GL32
10550/07	500	0.20	75	GL45
10550/08	500	0.45	75	GL45
10550/09	1000	0.20	90	GL32
10550/10	1000	0.20	90	GL45

PES Membrane (Blue collar)



Product no	Reservoir vol. (ml)	Pore size (µm)	Membrane Dia (mm)	Cap size
10551/01	500	0.20	75	GL32
10551/11	500	0.45	75	GL32
10551/02	500	0.20	75	GL45
10551/12	500	0.45	75	GL45

Bottle-Top Filter Units, Re-usable


These graduated, polysulphone tops are autoclavable for re-use, and incorporate a membrane support plate for 47mm or 50mm membranes. The polypropylene tubing adaptor can be stuffed with cotton for sterile venting. Suitable glass bottles (such as Duran brand) to fit the neck size should be purchased separately.



Product No	Volume (ml)	Neck Size
10554/01	250	GL32
10554/02	250	GL45
10554/03	500	GL32
10554/04	500	GL45

Filter Holders, Re-usable

Made from autoclavable polysulphone with graduated receiver. Tubing adaptor can be stuffed with cotton for sterile venting.



Product No	Reservoir Vol. (ml)	Receiver Vol. (ml)
10555/01	250	250
10555/02	500	500
10555/03	500	1000

Funnel Conversion Assembly

Converts Nalgene Filter Units into filter funnels. A funnel and vacuum gasket are included. The outer diameter of the gasket tapers from 40 to 28mm.

Product no 10557/01

Funnel Holders With Funnel

Made from polysulphone, these holders can be used with any filter flask or manifold that accepts a suitable rubber stopper with hole. An easy to use, no-twist vacuum gasket is included. (The outer diameter of the gasket tapers from 40 to 28mm). For use with 47 or 50mm membranes.



Product No	Capacity (ml)
10556/01	250
10556/02	500

Filter Holder Bottle Converters

An autoclavable polysulphone assembly that allows the conversion of reusable filter holders to bottle top filters, for use with Duran type glass media bottles.

Product No	Neck Size
10558/01	GL32
10558/02	GL45