

Flasks

Baffled Flasks

Borosilicate glass flasks with three deep baffles to promote turbulence.

Product no	Capacity (ml)
10888/01	100
10888/02	250
10888/03	500
10888/04	1000
10888/05	2000
10888/06	3000
10888/07	5000



Biological Stirrer Flasks

Borosilicate glass culture flasks each with three necks, the centre neck having a stirrer rod suspended from it. When used with a magnetic stirrer (see the "Stirrers" section) these culture flasks provide slow, gentle and even stirring of biological samples and cultures. Being borosilicate glass these flasks can be used in an incubator, and are supplied complete with heat resistant (up to 170°C) melamine caps and the stirring bar.



Product no	Total Volume (ml)	Working Volume	Centre Neck	Side Necks	Dia mm
10889/01	250	125	GL45	GL25	70
10889/02	500	250	GL45	GL25	86
10889/03	1000	500	GL45	GL25	101
10889/04	2000	1000	GL45	GL25	136
10889/05	5000	2500	GL45	GL45	186
10889/06	10000	5000	GL45	GL45	234
10889/07	20000	10000	GL45	GL45	300

Brewers Flasks

Glass Brewers Mash flasks to BS701 marked at 500 and 515 ml.

Product no
10594/01

Conical Flasks, Narrow Mouth

Erlenmeyer pattern flasks of heat-resistant borosilicate glass, with graduations, to ISO 1773, DIN 12380.

Product no	Capacity (ml)	Neck O.D. (mm)
10025/01	25	22
10025/02	50	22
10025/03	100	22
10025/04	200	34
10025/05	250	34
10025/06	300	34
10025/07	500	34
10025/08	1000	42
10025/09	2000	50
10025/10	3000	51
10025/11	5000	51



Conical Flasks, Heavy Duty

Designed for high mechanical strength with reinforced rims. Approximate graduations are shown in white enamel.

Product no	Capacity (ml)	Std Pack Quantity
11163/01	10	12
11163/02	25	48
11163/03	50	48
11163/04	125	48
11163/05	250	48
11163/06	500	36
11163/07	1000	24
11163/08	2000	8
11163/09	4000	4
11163/10	6000	4

Conical Flasks, Wide Mouth

Erlenmeyer pattern flasks of heat-resistant borosilicate glass, with graduations, to DIN 12385.



Product no	Capacity (ml)	Neck O.D. (mm)	Std Pack Quantity
10240/01	25	31	10
10240/02	50	34	10
10240/03	100	34	10
10240/04	200	50	10
10240/05	250	50	10
10240/06	300	50	10
10240/07	500	50	10
10240/08	1000	50	10
10240/09	2000	72	10