

Conical Flasks, Teflon® FEP

These translucent/clear Erlenmeyer flasks with Tefzel® screw caps are unbreakable with excellent chemical resistance. The flasks can withstand temperatures from -270°C to +150°C, while the caps can be used from -100°C to +150°C. Autoclavable. Supplied singly



Product no	Capacity (ml)
11164/01	50
11164/02	125
11164/03	250

Conical Flasks, Polycarbonate

Transparent polycarbonate Erlenmeyer flasks with polypropylene screw caps, offering excellent clarity. Autoclavable.



Product no	Capacity (ml)	Pack Size
11165/02	50	6
11165/03	125	6
11165/04	250	4
11165/05	500	4
11165/06	1000	2

Conical Flasks, PMP

Transparent polymethylpentene Erlenmeyer flasks with polypropylene screw caps, offering clarity, chemical resistance and impact strength. Autoclavable.



Product no	Capacity (ml)	Pack Size
11166/02	125	6
11166/03	250	4
11166/04	500	4
11166/05	1000	2

Culture Flasks


Flat-bottomed, thick walled borosilicate glass vessels for fermentation and culture work. Lids and Clamps can be found below and must be ordered separately.



Product no	Nominal Capacity (ml)	Flange Bore (mm)	Approx Height (mm)	Approx Dia (mm)
10684/01	500	75	125	85
10684/02	1000	100	160	105
10684/03	2000	100	225	130
10684/04	5000	100	290	185
10684/05	10000	100	380	215
10684/06	20000	100	500	250

Culture Flask Lids

Flat flange single and multi-socket glass covers for use with flasks (Product no 10684 above). Suitable flange clips can be found below (Product no 10683).



Product no	Flange Bore (mm)	Centre Socket	Parallel Socket	Side Sockets		
				5°	10°	15°
10682/01	75	19/26	14/23	19/26	19/26	19/26
10682/02	75	34/35	-	-	-	-
10682/03	100	19/26	19/26	2 x 19/26	19/26	-
10682/04	100	19/26	19/26	24/29	19/26	-
10682/05	100	19/26	14/23	14/23	14/23	29/32
10682/06	100	19/26	19/26	19/26	19/26	34/35
10682/07	100	24/29	19/26	19/26	19/26	34/35
10682/08	100	29/32	-	-	-	-
10682/09	100	29/32	14/23	14/23	14/23	29/32
10682/10	100	34/35	-	-	-	-
10682/11	100	34/35	19/26	-	2 x 19/26	-

Culture Flask Clamps

Simple wire clips for safely clamping ground glass flat flanges of the dimensions below.

Product no	Flange Bore (mm)
10683/01	75
10683/02	100

