

Flasks

Distillation Flasks

Made of heat-resistant borosilicate glass.

Product no	Capacity (ml)	Neck O.D. (mm)	Flask O.D. (mm)	Height (mm)	Side arm O.D. x length (mm)
10584/01	100	22	63	150	8x200
10584/02	250	34	85	200	8x200
10584/03	500	34	105	250	8x200
10584/04	1000	34	131	300	12x250



Distillation Flasks, Engler

borosilicate heat resistant glass distillation flasks of the capacities indicated below. The 125ml size conforms to ASTM D86. The side arms are angled at 75° to the vertical.

Product no	Capacity (ml)	Neck O.D. (mm)	Flask O.D. (mm)	Height (mm)	Side arm O.D. x length (mm)
10585/01	100	16	66	215	6 x 100
10585/02	125	17	68	215	7 x 100
10585/03	150	16	73	223	6 x 100



Filter Flasks

Erlenmeyer pattern vacuum flasks of heat-resistant borosilicate glass, with glass hose connection, to ISO 6556.


Product no	Capacity (ml)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)	Std Pack Quantity
10577/01	100	65	24	100	10
10577/02	250	85	35	155	10
10577/03	500	105	34	185	10
10577/04	1000	135	45	230	10
10577/05	2000	166	60	255	1



Similar to 10577 Range, but with plastic hose connection and plastic tubulature:

The polypropylene assembly with silicone seal (10583/01) is replaceable and provides safer connection.


Product no	Capacity (ml)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)	Std Pack Quantity
10579/01	100	65	24	105	10
10579/02	250	85	35	155	10
10579/03	500	105	34	185	10
10579/04	1000	135	45	230	10
10579/05	2000	166	60	255	1



As 10579 Range, but with plastic coating:

Protects against mechanical damage. Will hold broken glass together if flask should implode. Autoclavable to 135°C.

Product no	Capacity (ml)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)	Std Pack Quantity
10581/01	100	65	24	100	10
10581/02	250	85	35	155	10
10581/03	500	105	34	185	10
10581/04	1000	135	45	230	10
10581/05	2000	166	60	265	1



Spare for 10579 & 10581 Ranges

Product no	Description
10583/01	Replacement hose assembly