

Filter Flasks, Large Capacity

Large vacuum flasks of heat-resistant borosilicate glass, with glass hose connection.



Product no	Capacity (litres)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)
10578/01	3	170	60	295
10578/02	5	185	70	360
10578/03	10	240	70	420
10578/04	15	255	75	500
10578/05	20	290	75	535

Similar to 10578 Range, but with plastic hose connection and plastic tubulature:

The polypropylene assembly with silicone seal (10583/01) is replaceable and provides safer connection.

Product no	Capacity (litres)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)
10580/01	3	170	60	295
10580/02	5	185	70	360
10580/03	10	240	70	420
10580/04	15	255	75	500
10580/05	20	290	75	535

As 10580 Range, but plastic coated:

Protects against mechanical damage. Will hold broken glass together if flask should implode. Autoclavable to 135°C.

Product no	Capacity (litres)	Max. O.D. (mm)	Neck O.D. (mm)	Height (mm)
10582/01	3	170	60	295
10582/02	5	185	70	360
10582/03	10	240	70	420
10582/04	15	255	75	500
10582/05	20	290	75	535

Spares for 10580, 10581 and 10582 flasks:

Product no	Description
10583/01	Replacement hose assembly

Filter Flasks, Large Capacity

Moulded in one piece for strength, a barbed adapter accepts 1/4" to 3/8" I.D. tubing on the 500ml and 1000ml sizes and 3/8" on the 2000ml size. Autoclavable. Use with silicone stoppers (to be ordered separately).

Product no	Capacity (ml)
11167/01	500
11167/02	1000
11167/03	2000



Accessory Stoppers, Silicone

Product no	Stopper Size	To Suit
10219/48	29	11167/01
10219/50	33	11167/02
10219/54	41	11167/03

Boiling Flasks

Boiling flasks of heat-resistant borosilicate glass with beaded rim, to ISO 1773, DIN 12347.




Product no	Capacity (ml)	Neck O.D. (mm)	Height (mm)	Std Pack Qty
10362/01	50	26	95	10
10362/02	100	26	110	10
10362/03	250	34	144	10
10362/04	500	34	175	10
10362/05	1000	42	200	10
10362/06	2000	42	260	10
10362/07	3000	50	260	1
10362/08	4000	50	315	1
10362/09	5000	50	305	1
10362/10	6000	51	355	1
10362/11	10000	65	420	1
10362/12	20000	76	515	1

Flask Support

Made from autoclavable polypropylene, this support will safely hold flasks up to 4 litres capacity.

Product no
10838/01



Flask Rings, Cork

45mm depth composite compressed cork rings are highly heat and chemical resistant for use in supporting flasks.

Product no	Inside Dia (mm)	Outside Dia (mm)
10715/01	45	76
10715/02	85	112
10715/03	120	150
10715/04	160	195
10715/05	200	240