

Flasks

Flat Bottom Flasks, Narrow Neck

Made of heat-resistant borosilicate glass with beaded rim, to ISO 1773, DIN 12347.



Product no	Capacity (ml)	Neck O.D. (mm)	Height (mm)	Std Pack Qty
10407/01	50	22	100	10
10407/02	100	22	110	10
10407/03	250	35	140	10
10407/04	500	34	171	10
10407/05	1000	42	200	10
10407/06	2000	42	250	10
10407/07	3000	50	250	1
10407/08	4000	50	300	1
10407/09	5000	50	290	1
10407/10	6000	65	340	1
10407/11	10000	65	400	1

Kjeldahl Flasks

Made of heat-resistant borosilicate glass, to DIN 12360.



Product no	Capacity (ml)	Neck O.D. (mm)	Flask O.D. (mm)	Height (mm)
10576/01	50	22	51	200
10576/02	100	22	60	200
10576/03	250	34	81	270
10576/04	500	34	101	300
10576/05	750	34	115	340
10576/06	1000	34	126	350

Reaction Flasks

These flasks are the finest available, significantly better than conventional "flat flange" designs:

- Manufactured from heat-resistant borosilicate glass
- Schott flanges are used for their higher strength
- A grooved flange accepts an O-ring for a superior greaseless seal

Flange clamps:

- Allow instant assembly and disassembly
- Manufactured from stainless steel for strength and chemical resistance
- Incorporate a tension adjustment facility for a perfect fit

O-rings are available in standard silicone or FEP-coated rubber combining elasticity with excellent chemical resistance. Flange clamps and O-rings should be ordered separately.

Product no	Volume (ml)	Flange size (Bore-mm)
10692/01	1000ml	DN100
10692/02	2000ml	DN100
10692/03	3000ml	DN120
10692/04	4000ml	DN120
10692/05	5000ml	DN150
10692/06	6000ml	DN150
10692/07	10000ml	DN150
10692/08	20000ml	DN150

Accessories

Product no	Description	Flange Size	Pack Size
10177/01	Quick release clamp	DN60	1
10177/02	Quick release clamp	DN100	1
10177/03	Quick release clamp	DN120	1
10177/04	Quick release clamp	DN150	1
10177/05	Quick release clamp	DN200	1
10179/01	O-rings, silicone	DN60	5
10179/02	O-rings, silicone	DN100	5
10179/03	O-rings, silicone	DN120	5
10179/04	O-rings, silicone	DN150	5
10179/05	O-rings, silicone	DN200	5
10180/01	O-rings, FEP-coated	DN60	1
10180/02	O-rings, FEP-coated	DN100	1
10180/03	O-rings, FEP-coated	DN120	1
10180/04	O-rings, FEP-coated	DN150	1
10180/05	O-rings, FEP-coated	DN200	1