

Reaction Flasks, Flat Flange


Manufactured from borosilicate heat resistant glass. A wide access via the flat flange joint for ease of cleaning, loading and removal of solids etc, and for use with rigid paddle stirrers. Covers (listed below) are fully interchangeable and therefore one flask may suit many different reactions with just a change of lid.

Manufactured by Yorlab



Product no	Nominal Capacity (ml)	Flange Bore (mm)	Flange Shape	Approx Height (mm)
10681/01	250	75	Spherical	90
10681/02	500	75	Spherical	115
10681/03	700	100	Cylindrical	135
10681/04	1000	100	Spherical	155
10681/05	2000	100	Spherical	185
10681/06	3000	100	Spherical	210
10681/07	4000	100	Spherical	225
10681/08	5000	100	Spherical	245
10681/09	6000	100	Spherical	260
10681/10	10000	100	Spherical	320
10681/11	20000	100	Spherical	380

Flask lids



Product no	Flange Bore	Centre Socket	Side Sockets			
			Parallel	5	10	15
10682/01	75mm	19/26	14/23	19/26	19/26	19/26
10682/02	75mm	34/35	-	-	-	-
10682/03	100mm	19/26	19/26	2x19/26	19/26	-
10682/04	100mm	19/26	19/26	24/29	19/26	-
10682/05	100mm	19/26	14/23	14/23	14/23	29/32
10682/06	100mm	19/26	19/26	19/26	19/26	34/35
10682/07	100mm	24/29	19/26	19/26	19/26	34/35
10682/08	100mm	29/32	-	-	-	-
10682/09	100mm	29/32	14/23	14/23	14/23	29/32
10682/10	100mm	34/35	-	-	-	-
10682/11	100mm	34/35	19/26	-	2x19/26	-

Flange Clips (Wire):

Product no	Flange Bore (mm)
10683/01	75
10683/02	100

Flasks, Sugar Type 1

Two mark glass sugar flasks to BS675 in the capacities indicated below.

Product no	Graduation Marks
10592/01	100 & 110ml
10592/02	200 & 220ml



Flasks, Sugar Type 2/Kohlrausch

Glass sugar flasks suitable for sugar analysis. Other sizes available to order

Product no	Capacity (ml)
10593/01	100



Trypsinising Flasks

A range of sizes of trypsinising flasks all with three deep baffles to promote turbulence, and a dam below the sidearm to retain untrypsinised tissue when decanting. Neck sizes are standard throughout the range, and therefore the caps (which must be ordered separately) are fully interchangeable. Caps are anodized aluminium with an internally fitted chrome nickel spring which securely grips the inside of the tube providing a sterile but not hermetic seal.

Product no	Capacity (ml)	Neck O.D. (mm)	Side Arm O.D. (mm)
10887/01	250	38	13
10887/02	500	38	13
10887/03	1000	38	13
10887/04	3000	38	13
10887/05	5000	38	13

Accessory Caps:

Cap-O-Test Silver anodised aluminium caps with retaining spring. (Split packs available).

Product no	Size (mm)	Pack Size
10869/01	12/13	100
10869/51	37/39	50

Flask Brush

Curved bristle head 50mm dia x 100mm length. Joint between head and handle is hinged, and stem of handle is rubber coated to prevent breakages. Hardwood handgrip.

Product no
10125/01

