

# Flasks

## Volumetric Flasks, Class A

Manufactured to exacting standards and complying to BS 1792 - ISO 1042 (5000ml size is not BS/ISO listed). Each flask is numbered and individually calibrated. The flasks are manufactured from heat-resistant borosilicate glass EXCEPT the 15ml sizes (soda lime\*). A polypropylene stopper is supplied.

### Volumetric Flasks Without Works Certificate:

Product no	Capacity (ml)	Stopper Size	Std Pack Quantity
10221/01	5	10/13	5
10221/02	10	10/13	5
10221/03*	15	10/13	5
10221/04	20	10/13	5
10221/05	25	10/13	5
10221/06	50	10/13	5
10221/07	100	12/14	5
10221/08	200	14/15	2
10221/09	250	14/15	2
10221/10	500	19/17	2
10221/11	1000	24/20	2
10221/12	2000	29/32	1
10221/13	5000	34/35	1

\*soda lime glass



### Volumetric Flasks With Works Certificate:

Product no	Capacity (ml)	Stopper Size	Std Pack Quantity
10222/01	5	10/13	5
10222/02	10	10/13	5
10222/13*	15	10/13	5
10222/03	20	10/13	5
10222/04	25	10/13	5
10222/05	50	10/13	5
10222/06	100	12/14	5
10222/07	200	14/15	2
10222/08	250	14/15	2
10222/09	500	19/17	2
10222/10	1000	24/20	2
10222/11	2000	29/32	1
10222/12	5000	34/35	1

\*soda lime glass

### Amber Volumetric Flasks Without Works Certificate:

Product no	Capacity (ml)	Stopper Size	Std Packr Quantity
10223/01	5	10/13	5
10223/02	10	10/13	5
10223/03*	15	10/13	5
10223/04	20	10/13	5
10223/05	25	10/13	5
10223/06	50	10/13	5
10223/07	100	12/14	5
10223/08	200	14/15	2
10223/09	250	14/15	2
10223/10	500	19/17	2
10223/11	1000	24/20	2
10223/12	2000	29/32	1
10223/13	5000	34/35	1

\*soda lime glass

