

Amber Volumetric Flasks With Works Certificate:

Product no	Capacity (ml)	Stopper Size	Std Pack Quantity
11168/01	5	10/13	5
11168/02	10	10/13	5
11168/03*	15	10/13	5
11168/04	20	10/13	5
11168/05	25	10/13	5
11168/06	50	10/13	5
11168/07	100	12/14	5
11168/08	200	14/15	2
11168/09	250	14/15	2
11168/10	500	19/17	2
11168/11	1000	24/20	2
11168/12	2000	29/32	1
11168/13	5000	34/35	1

*Soda lime glass

Volumetric Flasks, Small Volume

Small volume Class A volumetric flasks in both clear and amber, with a trapezoidal shape for maximum stability. Works correction certificates are available on request.



Product no	Capacity (ml)	Accuracy (ml)	Colour
10961/01	1	±0.010	Clear
10961/02	2	±0.015	Clear
10961/03	3	±0.025	Clear
10961/04	4	±0.025	Clear
10961/05	5	±0.025	Clear
10961/11	1	±0.010	Amber
10961/12	2	±0.015	Amber
10961/13	3	±0.025	Amber
10961/14	4	±0.025	Amber
10961/15	5	±0.025	Amber

USP Class A Flasks

Each piece individually numbered and fully traceable satisfying USP/ISO/BS/DIN accuracies.

Class A - Clear

Product no	Capacity	USP	Pk Qty
11496/01	5ml	±0.02	5
11496/02	10ml	±0.02	5
11496/03	20ml	±0.03	5
11496/04	25ml	±0.03	5
11496/05	50ml	±0.05	5
11496/06	100ml	±0.08	5
11496/07	200ml	±0.10	2
11496/08	250ml	±0.12	2
11496/09	500ml	±0.15	2
11496/10	1000ml	±0.30	2

Volumetric Flasks, Plastic

Two ranges of plastic volumetric flask, both autoclavable. The polypropylene range will withstand autoclaving at 120°C whilst TPX will withstand 170°C and is crystal clear. Supplied singly with screw cap in the capacities indicated below.



Product no	Capacity (ml)	Material
10572/01	25	PP
10572/02	50	PP
10572/03	100	PP
10572/04	250	PP
10572/05	500	PP
10572/06	1000	PP
10572/11	50	TPX
10572/12	100	TPX
10572/13	250	TPX
10572/14	500	TPX
10572/15	1000	TPX