

Funnels

Glass Funnels, Economy

Simax borosilicate glass plain stem funnels. Supplied singly in the dimensions indicated below.

Product no	Top Dia (mm)
10288/01	55
10288/02	75
10288/03	100
10288/04	150
10288/05	200
10288/06	250



Borosilicate heat resistant glass short stem funnels. Supplied singly in the sizes indicated below.

Product no	Top Dia (mm)	Stem Length (mm)	Stem Dia (mm)	Filter Paper (mm)
10454/01	35	35	6	45-55
10454/02	45	45	6	55-70
10454/03	55	55	8	70-90
10454/04	70	70	8	110-125
10454/05	80	80	10	125-150
10454/06	100	100	10	150-185
10454/07	120	100	10	185-240
10454/08	150	150	16	240-270
10454/09	180	175	20	270-320
10454/10	200	175	26	320-385
10454/11	250	175	30	385-400
10454/12	300	175	30	500

Polypropylene Funnels

Each size takes a standard size filter paper and walls of body are angled at 60°. Ribbed inside and out to allow air to pass. Autoclavable to 120°C.



Product no	Top Dia (mm)	Stem Dia (mm)	Pack Size
11173/01	27	4	20
11173/02	37	5	20
11173/03	46	5	20
11173/04	66	10	20
11173/05	81	11	10
11173/06	100	11	1
11173/07	120	11	1
11173/08	150	14	1
11173/09	183	14	1

Polythene Funnels

Plain funnels in flexible polythene, with ribs for easier filling of liquids.

Product no	Dia (mm)
10374/01	37
10374/02	65
10374/03	90
10374/04	115
10374/05	140
10374/06	165
10374/07	205



H.D.P.E. Funnels

Tough, robust general purpose funnels. External ribs allow air to escape from receiver (not supplied).

Product no	Dia (mm)
10269/01	100
10269/02	120
10269/03	140
10269/04	180
10269/05	220
10269/06	260
10269/07	310



Buchner Funnels

Polypropylene two part Buchner funnels separate below the perforated disc for quick and easy cleaning. The robust construction makes these funnels virtually unbreakable and being made from polypropylene are resistant to most chemicals and temperatures to 120°C.



Product no	Capacity (ml)	Diaphragm Dia (mm)
10287/01	40	42.5
10287/02	70	55
10287/03	180	70
10287/04	285	80
10287/05	390	90
10287/06	810	110
10287/07	2100	160
10287/08	6000	240

Porcelain Buchner Funnels

Porcelain Buchner funnels with integral perforated support plate for the paper diameters below.

Product no	Paper Dia (mm)
10782/01	45
10782/02	55
10782/03	70
10782/04	90
10782/05	110
10782/06	125
10782/07	150
10782/08	185

Buchner Funnels Adaptor

A multi-purpose rubber ring that will adapt any size of funnel to a filter flask. **Product no 11169/01**