

Filter Funnels

Filter funnels with integral sintered glass disc made from heat resistant borosilicate glass. Supplied singly in the capacities and porosities indicated below. For porosity information see earlier pages. For larger funnels see Product no 10415 below.



Product no	Capacity (ml)	Disc Diam (mm)	Porosity	Stem O.D. (mm)
10414/01	50	35	1	10
10414/02	50	35	2	10
10414/03	50	35	3	10
10414/04	50	35	4	10
10414/05	50	35	5	10
10414/10	75	45	1	10
10414/11	75	45	2	10
10414/12	75	45	3	10
10414/13	75	45	4	10
10414/14	75	45	5	10
10414/20	125	60	1	10
10414/21	125	60	2	10
10414/22	125	60	3	10
10414/23	125	60	4	10
10414/24	125	60	5	10
10414/30	500	95	1	22
10414/31	500	95	2	22
10414/32	500	95	3	22
10414/33	500	95	4	22
10414/34	500	95	5	22
10414/40	1000	120	1	22
10414/41	1000	120	2	22
10414/42	1000	120	3	22
10414/43	1000	120	4	22
10414/44	1000	120	5	22

Filter Funnels, Larger Volumes

Similar to product no 10414 above but in larger sizes up to 30 litres and with a domed sintered disc which is more resistant to collapse under reduced pressure. Other sizes can be made to order.

Product no	Capacity (ml)	Disc Dia (mm)	Porosity	Stem O.D. (mm)
10415/01	3500	175	1	30
10415/02	3500	175	2	30
10415/03	3500	175	3	30
10415/04	3500	175	4	30
10415/10	5000	200	1	30
10415/11	5000	200	2	30
10415/12	5000	200	3	30
10415/20	7000	200	1	30
10415/21	7000	200	2	30
10415/22	7000	200	3	30
10415/30	10000	225	1	30
10415/31	10000	225	2	30
10415/32	10000	225	3	30
10415/40	12000	225	1	30
10415/41	12000	225	2	30
10415/42	12000	225	3	30
10415/50	20000	270	1	30
10415/51	20000	270	2	30
10415/52	20000	270	3	30
10415/60	30000	270	1	30
10415/61	30000	270	2	30
10415/62	30000	270	3	30

Nesting Stopper Guko Sets

Made from EDPM rubber, this set of eight gaskets will withstand up to 150°C.

Product no **Set of 8**
10840/01 22-84mm

Gasket	1	2	3	4	5	6	7	8
Large OD (mm)	22	29	36	44	53	63	73	84
Small OD (mm)	12	16	22	27	33	43	52	61
Small ID (mm)	7	9	15	19	24	33	42	50
Height (mm)	17	22	25	30	35	35	40	40

Powder Funnels

Polypropylene powder funnels have a wide bore short stem as indicated below for fast transfer of powders and liquids. All sizes have both internal and external ribs.

Product no	Dia (mm)	Stem Dia (mm)
10270/01	60	15
10270/02	80	15
10270/03	100	25
10270/04	120	30
10270/05	150	35
10270/06	180	40

Cream Weighing Funnels

Glass funnels with ground glass joint to BS696 for weighing cream into butyrometers.

Product no	Description
10573/01	Funnel only
10573/02	Stopper only

Funnel Holders

Single or double position polypropylene holders suitable for funnel diameters from 25mm to 170mm.

(Will fit retort rods from 8mm to 14mm.)

Product no	No of Positions
10833/01	1
10833/02	2

