

Dial Hygrometer

Direct reading dial instrument with reset pointer to indicate the last reading. Range of 0 to 100% humidity.

Product No	Dial colour	Movement type
10895/02	White	Spiral

Masons Pattern Hygrometer

Housed in a yellow polypropylene case measuring 380 x 80mm and complete with lid and bulb protectors and detachable water cistern. Range of -5°C to +50°C, (+20°F to +120°F). Includes hygrometric tables, spare wick and mounting screws.



Product No 10893/01

Whirling Hygrometer

Wet and dry whirling hygrometer, conforming to BS 2842. Range of -5°C to +50°C. Includes hygrometric tables and spare wick. A carry case with shoulder strap is available as an accessory.



Product No	Description
10894/01	Hygrometer
10894/02	Case

Hygrometer, Waterproof

This portable, heavy duty and water- and dust-proof meter takes Relative Humidity and temperature measurements (user selectable for °C or °F). A probe with two metre cable is supplied as standard or a probe with 5-metre cable is available for even greater versatility. The meter is supplied pre-calibrated so it is received ready for use.



RH range	5 to 95%
Temperature range	0 to +60°C, (+32 to 140°F)
Resolution	0.1% RH, 0.1°C, 1°F
Accuracy	±2%RH, ±0.4°C, ±1°F

Product No	Description
11228/01	Hygrometer
11228/11	RH probe, 5 metre cable
10993/01	Spare battery, 1.4 volt AA (Four batteries required)

Thermo-Hygrometer, Min/Max



Bench Mounting

A low-cost meter ideal for QC work, available with an optional ISO 9000 Calibration Certificate. A keyhole slot allows meter to be wall mounted or will stand upright on a bench. The meter gives simultaneous 23mm displays of humidity and temperature in separate LCD windows. Both temperature and humidity memories are continually updated and may be recalled at the touch of a button. (These values may then be reset by the user). A battery (size AAA) is required.

Product No	Description
10989/01	Thermo-hygrometer
10993/01	Battery, 1.5 volt AAA
10989/02	Calibration certificate (ISO 9000)

YORLAB