

Ice

Ice cube maker



Range of ice cube makers to suit the requirements of most laboratory applications. Interior and exterior stainless steel. Available in a range of throughputs, with or without storage facility. Electrical supply 230V, 50Hz.

Product no	Output, kg/24h	Storage, kg
12920/01	10	3.5
12920/02	24	14
12920/03	31	14
12920/04	39	14
12920/05	50	18
12920/06	71	28
12920/07	86	28
12920/08	145	70

**PLEASE ENQUIRE IF YOU
REQUIRE INSTALLATION**

Ice flake maker



Flake ice is stored inside and insulated plastic-lined storage bin with front opening for easy access. Electrical and water connections are located on the rear panel. The compressor is air cooled. Available in a range of sizes with or without storage facility. Interior and exterior made from stainless steel.

The above ice machines are intended for connection to a mains water supply. Connection to other sources may be possible, but may require extra parts. The installation of a point of use water filler is always recommended.

Standard services (within 1m or the machine's position) to be supplied by end use prior to delivery:

1. Correctly rated power socket
2. 3/4" washing machine tap on mains water supply
3. 1 1/2 " vented drain, no higher than base of washing machine

Note: 30mm required behind machine for utility connections
150mm space around machines for ventilation

Electrical supply 230V, 50Hz

Please enquire for details of installation and removal of old machines

Product no	Output, kg/24h	Storage, kg
12920/10	85	32
12920/11	70	25
12920/12	120	70
12920/13	200	70