

Receiver Adaptors



Short & Long Stem

Product No	Socket size	Approx stem length-mm
10645/01	14/23	65
10645/02	19/26	65
10645/03	24/29	65
10645/04	29/32	65
10646/01	10/19	120
10646/02	14/23	190
10646/03	19/26	200
10646/04	24/29	200
10646/05	29/32	200

Straight



Product No	Socket size	Overall height-mm
10648/01	19/26	180
10648/02	24/29	180

Vacuum, Straight

With a tubulatur for connection to a vacuum. Available with other joint sizes

Product No	Socket size	Cone size	Approx overall height-mm
10651/01	24/29	24/29	140



Plain Bend

Product No	Socket size	Cone size	Approx overall height-mm
10650/01	14/23	14/23	90
10650/02	19/26	19/26	100
10650/03	19/26	24/29	110
10650/04	24/29	24/29	125
10650/05	29/32	29/32	140
10650/06	34/35	34/35	155

Vented

Product No	Socket size	Cone size	Approx overall height-mm
10647/01	10/19	10/19	75
10647/02	14/23	14/23	80
10647/03	19/26	19/26	106
10647/04	19/26	24/29	120
10647/05	24/29	24/29	130
10647/06	29/32	29/32	140
10647/07	34/35	34/35	150

Vacuum, Angled

Angled receiver adaptor with a tubulatur, angled for connection to a vacuum supply.

Manufactured by Yorlab



Product No	Socket size	Cone size	Approx overall height-mm
10649/01	14/23	14/23	110
10649/02	19/26	19/26	140
10649/03	19/26	24/29	140
10649/04	24/29	24/29	150
10649/05	29/32	29/32	165
10649/06	34/35	34/35	180