


# Jointed Glassware

## Tubing Adaptors


### Right-angled, With Cone

Product No	Cone size	Tubing size-mm
10641/01	10/19	9
10641/02	14/23	9
10641/03	19/26	9
10641/04	24/29	9
10641/05	29/32	9
10641/06	34/35	9



### Straight, With Socket

Product No	Socket size	Tubing size-mm
10642/01	19/26	9
10642/02	24/29	9
10642/03	29/32	9



## Adaptors with Stopcock

Adaptors incorporating a GP Rotaflo stopcock and tubing connector. Provides gas or liquid feed or venting with fine control.



Product No	Cone size	OD of Olive-mm	Stopcock bore-mm
10644/01	14/23	9	3
10644/02	19/26	9	3
10644/03	24/29	9	3
10644/04	29/32	9	3

## Dreschel (Gas Washing) Heads

### Plain

Plain tube gas washing heads all with 24/29 cones. For use with bottles Cat No 10657.



Product No	For bottle capacity-ml	Approx stem length-mm
10658/01	125	130
10658/02	250	175
10658/03	500	215

### Adjustable:

Similar to Product No 10658 but with a screwthread joint on the receiving arm enables this bottle head to be adjusted for use with all sizes of bottle (Product No 10657).

Product No	Cone size	Thread size
10659/01	24/29	13

### Sintered

Similar to Cat No 10658 but available with a range of sinters at the end of the receiving arm to assist gas distribution. All with 24/29 cones.



Product No	Sinter porosity	For bottle capacity-ml	Stem length-mm
10660/01	0	125	132
10660/02	0	250	165
10660/03	0	500	217
10660/04	1	125	132
10660/05	1	250	165
10660/06	1	500	215
10660/07	2	250	165
10660/08	2	500	217
10660/09	3	250	165

### Bottles

For Dreschel bottles see section under Glass jars and bottles