


Jacketed Coil

Borosilicate glass coil condensers with a cooled outer surface. Supplied complete with screwthread tubing connectors.



Product No	Approx surface area x10-2m ²	Socket size	Cone size	Approx.eff ect. length-mm	Approx overall length-mm
10667/01	4.0	19/26	19/26	160	360
10667/02	4.0	24/29	24/29	160	365
10667/03	4.0	29/32	29/32	160	370

Liebig

Manufactured from borosilicate heat resistant glass, these condensers are water cooled and supplied with detachable screwthread connectors for safe and easy attachment of tubing (Cat No 10662/01 has glass olives only).



Product No	Approx surface area x10-2m ²	Socket size	Cone size	Approx. effect. length-mm	Approx overall length-mm
10662/01	0.15	10/19	10/19	60	135
10662/02	0.6	14/23	14/23	155	260
10662/03	1.0	19/26	19/26	208	340
10662/04	1.0	19/26	19/26	250	380
10662/05	1.0	24/29	24/29	250	390
10662/06	1.3	29/32	29/32	250	400
10662/07	2.0	19/26	19/26	380	530
10662/08	2.0	24/29	24/29	395	535
10662/09	2.9	29/32	29/32	400	550

Dean And Stark Apparatus

Dean and Stark receivers for determination of water content in light liquids. Manufactured to BS756 in the dimensions below. Cat No 10693/03 has a stopcock.



Product No	Capacity-ml	Graduations-ml	Socket size	Cone size
10693/01	2.0	0.05	19/26	24/29
10693/02	7.5	0.10	19/26	24/29
10693/03	10.0	0.10	19/26	24/29
10693/04	25.0	0.20	19/26	24/29

Soxhlet Extractors

Plain Body

Plain extractor bodies can be used with liners (product No **10673** below) or thimbles found in the Extraction section.



Product No	Nominal capacity-ml	Socket size	Cone size	Approx overall height
10672/01	60	34/35	24/29	205
10672/02	100	45/40	24/29	225
10672/03	200	55/44	24/29	250

Soxhlet Liners

Borosilicate glass syphon assembly which used to form a Jacketed Soxhlet.

Product No	Nominal capacity-ml *	To fit body
10673/01	20	10672/01
10673/02	50	10672/02
10673/03	100	10672/03

* Nominal capacity (siphoning volume)