

Jointed Glassware

BS2071

Manufactured from borosilicate heat resistant glass these Soxhlets are suitable for use with a wide range of solvents when used with paper thimbles, which can be found in the Extraction section.



YORLAB HAVE A CUSTOM GLASS BLOWING FACILITY. WE ARE ABLE TO MANUFACTURE TO YOUR OWN DESIGN

AN EXAMPLE OF OUR BESPOKE WORK BELOW IS CHITTICK APPARATUS FOR THE DETERMINATION OF CARBON DIOXIDE. THIS HAS BEEN MADE FOR SEVERAL OF OUR CUSTOMERS EACH TO THEIR OWN BESPOKE NEEDS.

Product No	Nominal capacity siphoning volume-ml	Socket size	Cone size	Approx.ove rall height-mm	Recon'd thimble size-mm
10670/01	6	19/26	14/23	150	10x50
10670/02	20	24/29	19/26	175	19x65
10670/03	40	29/32	24/29	200	22x80
10670/04	60	34/35	24/29	205	28x80
10670/05	60	34/35	34/35	205	28x80
10670/06	100	34/35	34/35	240	28x120
10670/07	100	40/38	24/29	220	33x100
10670/08	200	50/42	24/29	250	41x123

Large Capacity

Similar to Cat No 10670 but with a flat flange joint at the top of the body for easy loading. Suitable lids (10682 range) can be found on the following pages

Product No	Nominal capacity siphoning volume-ml	Flange size-mm	Cone size
10671/01	600	75	34/35
10671/02	2000	100	34/35

