

Reaction Flasks, Flat Flange



Manufactured from borosilicate heat resistant glass with wide access via the flat flange joint for ease of cleaning, loading and removal of solids etc, and use with rigid paddle stirrers. Covers (listed below) are fully interchangeable and therefore one flask may suit many different reactions with just a change of lid.

Product No	Nominal capacity-ml	Flange bore-ml	Flange shape	Approx height-mm
10681/01	250	75	Spherical	90
10681/02	500	75	Spherical	115
10681/03	700	100	Cylindrical	135
10681/04	1000	100	Spherical	155
10681/05	2000	100	Spherical	185
10681/06	3000	100	Spherical	210
10681/07	4000	100	Spherical	225
10681/08	5000	100	Spherical	245
10681/09	6000	100	Spherical	260
10681/10	10000	100	Spherical	320
10681/11	20000	100	Spherical	380

Flask Covers

Product No	Flange bore	Centre socket	Side sockets			
			Parallel	5"	10"	15"
10682/01	75	19/26	14/23	19/26	19/26	19/26
10682/02	75	34/35	--	--	--	--
10682/03	100	19/26	19/26	2x19/26	19/26	--
10682/04	100	19/26	19/26	24/29	19/26	--
10682/05	100	19/26	14/23	14/23	14/23	29/32
10682/06	100	19/26	19/26	19/26	19/26	34/35
10682/07	100	24/29	19/26	19/26	19/26	34/35
10682/08	100	29/32	--	--	--	--
10682/09	100	29/32	14/23	14/23	14/23	29/32
10682/10	100	34/35	--	--	--	--
10682/11	100	34/35	19/26	--	2x19/26	--

Flange Clips (Wire):

Product No	Quickfit® alternative	Flange bore-mm
10683/01	JC75F	75
10683/02	JC100F	100

