

# Measures and Microbiology

## Imhoff Cone

Made of clear, autoclavable polycarbonate with polypropylene leak-proof screw closure. Fully graduated with divisions:

0-1ml x 0.1ml,  
1-10ml x 0.5ml,  
10-20ml x 1ml,  
20-40ml x 2ml,  
40-100ml x 10ml,  
100-500ml x 50ml and  
500-1000ml x 100ml.  
Top diameter is 106mm.



Product No	Description
10908/01	Imhoff cone
10909/01	Stand - 2 position
10909/02	Stand - 3 position

## Jugs

### Polypropylene

Graduated, jugs, highly translucent and autoclavable.

Product No	Capacity-ml	Grads-ml
10461/01	500	25
10461/02	1000	50
10461/03	2000	50
10461/04	3000	100
10461/05	5000	250



### HDPE

Beaker jugs, graduated in millilitres and ounces. Made of heavy-wall HDPE, and suitable for acids, alkalis, etc.


Product No	Capacity-ml
11244/01	1000
11244/02	3000



## Dry Swab

Product No 11256/11

Blue plasticised 12" paper shaft with jumbo dry foam tip, in flat-bottom tube. Gamma irradiated. Supplied in packs of 100.



## Amies Transport Swabs

Gamma irradiated Probact swabs in a round bottom tube of transport medium in a peel pouch. A hexagonal cap prevents rolling. Supplied in packs of 100

### Amies Clear in Peel Pouch

Product No	Description
10851/01	Woodstem, cotton tip
10851/10	Aluminium shaft, cotton mini tip
10851/02	Polystyrene shaft, pure viscose tip

### Amies Charcoal in Peel Pouch

Product No	Description
10851/03	Woodstem, cotton tip
10851/04	Polystyrene shaft, pre viscose tip
10851/05	Fine twisted wire stem, cotton mini tip

## Cotton Wool Non-Sterile

Supplied in 500g rolls in a range of colours and absorbent or non-absorbent as shown below. Ideal for plugging culture tubes for a breathing plug

Product No	Description
10233/01	Non-absorbent white
10233/02	Non-absorbent blue
10233/03	Non-absorbent green
10233/04	Non-absorbent yellow
10233/05	Non-absorbent red
10233/06	BP quality, absorbent
10233/07	Hospital quality, absorbent

### Cotton Wool Balls

Product No	Description
10048/01	Standard quality, pack of 100