


# Microbiology

## Petri Dishes

### Sterile Polystyrene


Aseptically manufactured in crystal clear polystyrene. Contract prices available for large users.



| Product No | Diameter-mm | No. of vents/style | Case qty | Sleeve qty |
|------------|-------------|--------------------|----------|------------|
| 10268/03   | 55          | 0                  | 1620     | 15         |
| 10268/04   | 55          | 3                  | 1620     | 15         |
| 10268/05   | 55          | Contact plate      | 720      | 15         |
| 10268/06   | 90          | 0                  | 600      | 25         |
| 10268/01   | 90          | 1                  | 500      | 20         |
| 10268/02   | 90          | 3                  | 500      | 20         |
| 10268/07   | 90          | 2 compartment      | 500      | 20         |
| 10268/08   | 140         | 3                  | 176      | 16         |
| 10268/60   | 90          | 1                  | 500      | 25         |
| 10268/61   | 90          | 3                  | 500      | 25         |

### Glass

Manufactured to DIN 12339 and allow even distribution of media. The borosilicate glass base and lid have distortion-free transparency.



| Product No | Base od x height*-mm |
|------------|----------------------|
| 10608/01   | 54.5 x 22            |
| 10608/02   | 74 x 22              |
| 10608/03   | 94 x 17              |
| 10608/04   | 94 x 22              |
| 10608/05   | 114 x 22             |
| 10608/06   | 143 x 32             |

\*height includes lid

## Culture/Roller Bottles

Made of borosilicate glass. Other sizes available. Caps (GL45) should be ordered separately

| Product No | Outside dia-mm | Length-mm | Wall thickness-mm |
|------------|----------------|-----------|-------------------|
| 10890/01   | 110            | 350       | 5                 |
| 10890/02   | 110            | 550       | 5                 |

### Accessory Caps

GL45 caps autoclavable to 140°C. Supplied in packs of 10.

| Product No | Colour |
|------------|--------|
| 10020/03   | Blue   |
| 10020/04   | Yellow |
| 10020/05   | Green  |
| 10020/06   | Grey   |
| 10020/07   | White  |

### Duran GL45 Red Caps with PTFE Seal

Will withstand dry heat sterilisation up to 200°C. Supplied in pack of 10

| Product No | Colour |
|------------|--------|
| 10021/03   | Red    |

## Culture Loops

Sterile disposable, supplied in boxes of 1000, inner sleeves of 20.

| Product No | Volume (µl) |
|------------|-------------|
| 10030/01   | 1           |
| 10030/03   | 10          |