

Microscopes

XEC 80 and 85 Series

Professional Stereo Microscopes




XES-80-2L Zoom Stereo Microscope

The XES series is a range of professional stereo microscopes for use in both general and industrial fields. Interchangeable or zoom lenses, together with high quality optics ensure good image quality throughout the range.

The XES-80 stereo microscopes have a 1:4 zoom range and flexibility to suit all types of application. Both binocular and trinocular head versions are available and all models feature independently adjustable transmitted and incident lighting. The XES-85 is a modular series of stereo zoom microscopes and systems.



XES-80 consists of the following:

Body	Die cast aluminium
Head	45o inclined binocular head
Eyepiece	Paired WF 10x/20 with ± 5 diopter adjustment
Objective	Zoom range 1:4 (7.5x..35x..70x with optional objectives)
Stage/Stand	110mm x 120mm with sprung clips
Focusing	Rack and pinion
Illumination	Halogen Incident and and transmitted lighting



Order Information

Product No	Description
13350/20	XES-80-2L 1:4 Stereo Zoom Microscope
13350/21	XES-80T-2L 1:4 Stereo Zoom Microscope with trinocular head and CCTV adaptor.
13350/22	Binocular 1:4 Stereo Zoom Head
13350/23	Binocular 1:6 Stereo Zoom Head
13350/24	Trinocular 1:4 Stereo Zoom Head
13350/25	Plain stand with rack and pinion focusing mount
13350/26	5-axis Stand without focusing mount with 150mm vertical post and 200mm horizontal post
13350/27	Plain stand without focusing mount with 150mm vertical post and 200mm horizontal post
13350/28	Rack and pinion focusing mount for 723-107/108



Binocular head, rack and pinion focusing mount, plain stand with vertical and horizontal posts.

Many other accessories available