

Microscope Accessories

Lens & Optical Glass Cleaning

We offer a selection of lens tissues and wipes:



Product No	Brand & Grade	Description	Size-mm	Pack size
12350/01	Schleicher & Schuell	Thin, white & non-fibrous Tissue paper	500 x 760	500 sheets
12350/02	2478		100 x 150 (special shapes)	100
10162/01	Kleenex 3020	Low-lint tissues in a splash proof box	210 x 110	144 boxes of 76 tissues
10175/01	Whatman 105	Highly absorbent lens tissues with a soft surface	100 x 150	25 wallets of 25 sheets
10175/02			200 x 300	100 sheets
10175/03			460 x 570	500 sheets


Coverslips

Glass coverslips conforming to BS 7011. Thickness No. 1 is 0.13-0.17mm and No. 1.5 is 0.16-0.19mm. Supplied in plastic boxes of 100.

Product No	Thickness	Size-mm	Product No	Thickness	Size-mm
10747/01	No.1	18 x 18	10747/11	No.1.5	18 x 18
10747/02	No.1	22 x 22	10747/12	No.1.5	22 x 22
10747/03	No.1	22 x 26	10747/13	No.1.5	22 x 26
10747/04	No.1	22 x 32	10747/14	No.1.5	22 x 32
10747/05	No.1	22 x 40	10747/15	No.1.5	22 x 40
10747/06	No.1	22 x 50	10747/16	No.1.5	22 x 50
10747/07	No.1	22 x 64	10747/17	No.1.5	22 x 64
10747/08	No.1	24 x 32	10747/18	No.1.5	24 x 32
10747/09	No.1	24 x 40	10747/19	No.1.5	24 x 40
10747/10	No.1	24 x 50	10747/20	No.1.5	24 x 50

Microscope Slides

Conforming to BS 3836, these slides are of soda glass, measure 76x26mm and have a nominal thickness of 1mm. The Twin Frost versions have writing areas on both sides. Supplied in packs of 50.



Product No	Description
10746/01	Plain unground
10746/02	Twin frost, unground
10746/03	Plain, ground 45°
10746/04	Plain, ground 90°
10746/05	Twin frost, 45°ground
10746/06	Twin frost, 90°ground

Cavity Slides

These pre-washed slides are of extra white glass, measuring 76x26mm with a nominal thickness of 1.5mm. The edges are ground. Cavities measure 15mm diameter and have a depth of 0.6-0.8mm. The 2-cavity slide is suitable for melting point determinations. Supplied in boxes of 50.

Product No	No. of cavities
10748/01	1
10748/02	2
10748/03	3