

Microscope Accessories

Fibre Optic Lamps



Provide light just where you need it, without heat. Select a self-supporting light guide most suited to the application.

- Vibration free
- Compact
- Maintenance-free
- Variable levels of brightness with mechanical dimming facility
- Bulb changing is easy and safe
- Stackable to save space when more than one lamp used
- KL 2500 has integral colour wheel and remote control option

Technical data

| Product No | 11267/01 | 11268/01 | 11268/02 |
|------------------------------|-------------------------|---------------|---------------|
| | KL200 | KL1500 LCD | KL2500 LCD |
| Dimensions (H x W x D-mm) | 120x125x205 | 170x200x265 | 170x200x265 |
| Approx weight-kg | 1.5 | 4.8 | 6.0 |
| Cooling | Natural convection | Low noise fan | Low noise fan |
| Power supply * | All 230 volts, 50/60 Hz | | |

Bulb

| Type | All halogen reflector | | |
|-------------------|-----------------------|-----------|-----------|
| Voltage | 8 | 15 | 24 |
| Power consumption | 20 watts | 150 watts | 250 watts |

Average life

| | | | |
|-------------------|-----|------|------|
| Stage 3 - (hours) | 300 | -- | -- |
| Stage 4 - (hours) | -- | 1500 | 1500 |
| Stage 5 - (hours) | -- | 150 | 150 |

* Power supply : Models to suit other supplies available

Natural Convection Cooling- KL200

| Product No | Model |
|------------|--------------------------------------|
| 11267/01 | KL200 light source |
| 11267/11 | Goose neck light guide, 1 arm, 500mm |
| 11267/12 | Goose neck light guide, 2 arm, 500mm |
| 11267/13 | 6-point ringlight, 58mm, 900mm |
| 11267/14 | Spare bulb |

Fan Cooling- KL1500 LCD & KL2500 LCD

| Product No | Model |
|------------|--------------------------------------|
| 11268/01 | KL 1500 LCD light source |
| 11268/02 | KL 2500 LCD light source |
| 11268/11 | Goose neck light guide, 1 arm, 600mm |
| 11268/12 | Goose neck light guide, 2 arm, 600mm |
| 11268/13 | Goose neck light guide, 3 arm, 600mm |
| 11268/14 | Continuous ringlight, 58mm, 600mm |