

# Mixing Heating and Shaking

## Stirrers, overhead EUROSTAR cont.

### Technical specifications

	<b>13920/20</b>	<b>13920/21</b>	<b>13920/22</b>	<b>13920/23</b>	<b>13920/24</b>	<b>13920/25</b>
	<b>Digital</b>	<b>Power basic</b>	<b>Power control visc</b>	<b>Power control visc, P1</b>	<b>Power control visc, P4</b>	<b>Power control visc, P7</b>
Stirring capacity, L (H <sub>2</sub> O)	20	40	40	60	40	40
Output at stirring shaft, W	53	105	105	126	95	95
Speed range, rpm	50 to 2,000	50 to 2,000	50 to 2,000	50 to 1,200	14 to 530	8 to 290
Speed display	Digital	Analogue	Digital	Digital	Digital	Digital
Dimensions (w x d x h)mm	80 x 190 x 222	80 x 190 x 253	80 x 190 x 205	80 x 90 x 253	80 x 190 x 253	80 x 190 x 330
Mass, kg	2.8	3.8	3.8	4	4.9	4.9
Electrical supply	230V	230V	230V	230V	230V	230V

### Accessories and spares

#### Product no Description

<b>13920/36</b>	Plate stand, 800mm height
<b>13920/37</b>	Boss head clamp
<b>13920/38</b>	Turbine stirrer, 50mm length x 8mm dia.
<b>13920/39</b>	Strap clamp
<b>13920/40</b>	Shaft protection, Plexiglas, length and adjustment 190/310mm

### Accessories EUROSTAR power basic and control visc.

#### Product no Description

<b>13920/50</b>	Telescopic stand, 620mm to 1,010mm
<b>13920/51</b>	Boss head clamp
<b>13920/52</b>	Paddle stirrer 70mm length x 8mm dia
<b>13920/53</b>	Anchor stirrer
<b>13920/54</b>	Strap clamp
<b>13920/55</b>	Shaft protection, Plexiglas, length adjustment 190/310mm

### Extra accessories EUROSTAR power control visc.

#### Product no Description

<b>13920/56</b>	PC 1.5m cable, 2.5m length
<b>13920/57</b>	Labworldsoft software
<b>13920/58</b>	AM 1 analogue adapter