

Mixing, Heating and Shaking

Stirrer, RW 28 basic



Powerful, mechanically controlled stirrer. Suitable for quantities up to 80 l.

- **Two speed ranges for highly viscous media and for intensive mixing**
- **Push-through agitator elements**

Laboratory stirrer for simple stirring tasks of up to 10L (H₂O) with ideal speed range from 40 to 1,200rpm. Especially suitable for schools, universities and inspection laboratories.

Technical specifications

Product No 13920/59

Motor rating input/output, W	90/220
Output at stirring shaft, W	90
Speed range, rpm	60-1400
Chuck range, mm	1-10
Diameter/length of extension arm, mm	145
Dimensions (w x d x h), mm	123 x 364 x 252
Mass, kg	7.4
Electrical supply, V	220-240

Product no Description

13920/60	Plate stand, 800mm height
13920/61	Boss head clamp
13920/62	Strap clamp
13920/63	Shaft protection, Plexiglas, length and adjustment 190/310mm
13920/64	Paddle stirrer 70mm length x 8mm dia
13920/65	Anchor stirrer, 350mm length x 8mm dia.
13920/66	4 bladed propeller stirrer, 350mm length x 8mm dia
13920/67	Shaft protection, Plexiglas, length adjustment 190/310mm

RW 20 Digital



- **With digital display**

- **robust, slimline, ergonomic design**

- **Technical improvements on the trusted RW 20 series designs**

- **With constant power drive**

- **Two speed ranges for universal use from 60 - 2.000 rpm**

- **Push-through agitator shafts (only when stationary)**

Technical specifications

Product No 13920/71

Speed range I (50Hz), rpm	RW20.n 60 to 500
Speed range II (50Hz), rpm	240 to 2,000
Motor rating input/output, W	70/35
Output stirring shaft	26
Chuck range, mm	0.5 to 10
Dimensions (w x d x h), mm	88 x 294 x 212
Mass, kg	3.1
Electrical supply, V	230-240

Product no Description

13920/72	Stirrer motor
13920/73	Digital

Accessories and spares

13920/80	Plate stand, 800mm height
13920/81	Boss head clamp
13920/82	Strap clamp
13920/83	Shaft protection, Plexiglas, length and adjustment 190/310mm
13920/84	Paddle stirrer 70mm length x 8mm dia
13920/85	Anchor stirrer, 350mm length x 8mm dia.
13920/86	4 bladed propeller stirrer, 350mm length x 8mm dia