

Mixing, Heating and Shaking

IKA Horizontal shaker HS 260 basic

Compact, flat horizontal shaker with ideal swivel motion.

- Electronic adjustment of speed and timer.
- LED display for speed and time adjustment.
- Attachments are not included in delivery

Accessories

-Attachments AS 260.1, AS 260.2, AS 260.3, AS 260.5

Shaking motion:	reciprocating
Stroke	20mm
Max. shaking weight: (incl. attachment)	7.5kg
Motor rating input/output:	45/10W
Adjustable Speed range:	20-300 min⁻¹
Speed display:	LED line
Timer (digital):	∞/5 - 50min
Dimensions (w x d x h)	360 x 420 x 100mm
Weight:	8.5kg
Voltage:	230V. 50/60Hz
	115V - on request!
Interface:	RS 232/analog
Protection class:	IP 21



Product No	Description
13960/01	HS 260 basic

IKA Horizontal shaker HS 501digital

Low profile horizontal laboratory shaker with a pleasant design and large mounting surface.

- Infinitely variable speed control
- Ideal for all lying vessels, e.g. separating funnels
- Guaranteed continuous operation even under extreme loads
- Includes timer
- Attachments are not included in delivery

Accessories

-Attachments AS501.1, As 501.2, As 501.4, As 501.5

Shaking motion:	reciprocating
Stroke	30mm
Max. shaking weight: (incl. attachment)	15kg
Motor rating input/output:	70/19W
Adjustable Speed range:	0-300 min⁻¹
Speed display:	Digital
Timer (digital):	1/56min
Dimensions (w x d x h)	505 x 585 x 120mm
Weight:	26 kg
Voltage:	230V. 50/60Hz
	115V - on request!
Protection class:	IP 21



Product No	Description
13960/02	HS 501 digital