

Mixing Heating and Shaking

Yorlab are the
Official agents for Service and
repair of all IKA units

IKA Orbital shaker KS 130 control

Quiet orbital shaker ensures long life with ideal swivel motion (wear-free).

- Electronic adjustment of speed and timer
- Digital display for speed, timer function and operating modes.
- Electronic timer, also for continuous operation
- With integrated end point positioning (for automated robot-controlled sampling)
- all functions can be controlled and documented with software labworldsoft®
- Attachments are not included in delivery.

Accessories

- Attachments AS 130.1, AS 130.2, AS 130.3, AS 130.4
- labworldsoft®: PC 1.5 cable

Shaking motion:	orbital
Stroke	4mm
Max. shaking weight: (incl. attachment)	2kg
Motor rating input/output:	45/10W
Adjustable Speed range:	0-800 min⁻¹
Speed display:	Digital
Timer (digital):	1 - 599min
Dimensions (w x d x h)	270 x 316 x 98mm
Weight:	9.8kg
Interface	RS 232/analog
Voltage:	230V. 50/60Hz 115V - on request!
Protection class:	IP 21

Product No	Description
13960/09	KS 130 control

