

pH Testers

pH meters, bench, models 3505, 3510 and 3520



Technical Specification

		13620/20	13620/21 and 13620/22
Range	pH	-2 to +16	-2.000 to +19.999
	mV	-1.999 to +1,999	-1.999 to +1,999
	°C	-10 to +105	-10 to +105
Resolution	pH	0.01	0.001/0.01
	mV	1	0.1
	°C	0.1	0.1
Accuracy	pH	±0.02	±0.003
	mV	±1	±0.2
	°C	±0.5	±0.5
Dimensions (w x d x h mm)		210 x 250 x 60	210 x 250 x 60
	Electrical supply	9V d.c. ±10% at 230V, 50/60Hz	



Product no	Description
13620/20	1
13620/21	2
13620/22	3

Accessories

Product no	Description
13620/23	Combination electrode, epoxy
13620/24	pH electrode, glass bodied, combination
13620/25	Power supply
13620/26	Dust cover
13620/27	Interface cable kit for 13620/21 and 13620/22

