

# Pipettes

VOLAC graduated and bulb one-mark pipettes have markings in permanent amber stain, colour codes are fused ceramic and tips are carefully formed, ground and bevelled to ensure proper flow rates, and "TUFFTIP" treated for greater resistance to chipping and breakage. Sizes 5ml and above are constricted for cotton plugging. Class A pipettes are manufactured in accordance with BS 700 and ISO 835. They are individually calibrated, tested, numbered and timed, and are fully traceable. VOLAC pipettes comply with BS 1583 and ISO 648 for accuracy and flow rate.



## Bulb Pipettes, Class A

NC denotes verification only. WAC denotes pipettes with a Works Accuracy Certificate. All are individually calibrated, tested, marked with the time of flow to deliver and numbered for full traceability.

### Technical Data - WITHOUT Certificate

Product No	Capacity-ml	Colour code	Accuracy Limits-ml	Flow time to delivery-seconds	Std. pack quantity
10225/01	1	Blue	±0.008	10 - 20	5
10225/02	2	Orange	±0.01	10 - 25	5
10225/03	3	Black	±0.015	10 - 25	5
10225/04	4	Red	±0.015	10 - 25	5
10225/05	5	White	±0.015	15 - 30	5
10225/06	10	Red	±0.02	15 - 40	5
10225/07	15	Green	±0.025	15 - 40	5
10225/08	20	Yellow	±0.03	25 - 50	5
10225/09	25	Blue	±0.03	25 - 50	5
10225/10	50	Red	±0.05	30 - 60	5
10225/11	100	Yellow	±0.08	40 - 60	2

### Technical Data - WITH Certificate:

Product No	Capacity-ml	Colour code	Accuracy Limits-ml	Flow time to delivery-seconds	Std. pack quantity
10226/01	1	Blue	±0.008	10 - 20	5
10226/02	2	Orange	±0.01	10 - 25	5
10226/03	3	Black	±0.015	10 - 25	5
10226/04	4	Red	±0.015	10 - 25	5
10226/05	5	White	±0.015	15 - 30	5
10226/06	10	Red	±0.02	15 - 40	5
10226/07	15	Green	±0.025	15 - 40	5
10226/08	20	Yellow	±0.03	25 - 50	5
10226/09	25	Blue	±0.03	25 - 50	5
10226/10	50	Red	±0.05	30 - 60	5
10226/11	100	Yellow	±0.08	40 - 60	2

## Bulb Pipettes, Class B

Colour coded one mark bulb pipettes manufactured to BS1583. Permanent amber markings, and Volac's unique tempered "Tufftip" for durability.

Product No	Capacity-ml	Colour code	Accuracy Limits-ml	Flow time to delivery-seconds	Std. pack quantity
10224/01	1	Blue	±0.015	5 - 20	5
10224/02	2	Orange	±0.02	5 - 25	5
10224/03	3	Black	±0.03	5 - 25	5
10224/04	4	Red	±0.03	5 - 25	5
10224/05	5	White	±0.03	7 - 30	5
10224/06	10	Red	±0.04	8 - 40	5
10224/07	15	Green	±0.05	8 - 40	5
10224/08	20	Yellow	±0.06	9 - 50	5
10224/09	25	Blue	±0.06	10 - 50	5
10224/10	50	Red	±0.1	13 - 60	5
10224/11	100	Yellow	±0.15	26 - 60	2

### Pipettes, Red Label

These Class B pipettes offer excellent accuracy and outstanding value for money. They are easily distinguished by their permanent red enamel markings.

Product No	Capacity-ml	Std. pack quantity
10941/01	1	5
10941/02	2	5
10941/03	5	5
10941/04	10	5
10941/05	25	5
10941/06	50	5
10941/07	100	2
10941/08	20	5