

## Bulb Pipettes, Polypropylene

Unbreakable Class B polypropylene one mark bulb pipettes. Autoclavable and particularly suitable for use in food



Product No	Capacity-ml
10678/01	1
10678/02	2
10678/03	5
10678/04	10
10678/05	25
10678/06	50



## Graduated Pipettes, Type 1

Class A and Class B pipettes, calibrated zero at top, total capacity marked at shoulder. Non blow-out glass.



Product No	Class	Capacity-ml	Sub Div-ml	Colour code
10232/01	A	1	0.01	Yellow
10232/02	A	2	0.02	Black
10232/03	A	5	0.05	Red
10232/04	A	10	0.1	Orange
10232/05	A	25	0.2	Green
10227/01	B	1	0.01	Yellow
10227/02	B	2	0.02	Black
10227/03	B	5	0.05	Red
10227/04	B	10	0.1	Orange
10227/05	B	25	0.2	Green

Data (continued)

Product No	Accuracy Limits-ml	Flow time to delivery-seconds
10232/01	±0.006	7 - 10
10232/02	±0.01	8 - 12
10232/03	±0.03	10 - 14
10232/04	±0.05	13 - 17
10232/05	±0.10	15 - 21
10227/01	±0.01	2 - 10
10227/02	±0.02	2 - 12
10227/03	±0.05	5 - 14
10227/04	±0.10	5 - 17
10227/05	±0.20	9 - 21