

Pipettes

Serological Pipettes, Disposable Glass

Non-sterile Packs

Product No	Jet	Capacity	Grads	Length	Inner sleeve	Pack Qty
11316/01	Short	1ml	0.01ml	230mm	10	500ml
11316/02	Short	1ml	0.1ml	230mm	10	500ml
11316/03	Drawn	1ml	0.01ml	270mm	10	500ml
11316/04	Short	2ml	0.02ml	260mm	10	500ml
11316/05	Drawn	2ml	0.02ml	270mm	10	500ml
11316/06	Drawn	5ml	0.05ml	300-310mm	10	500ml
11316/07	Drawn	5ml	0.1ml	Shorty	5	250ml
11316/08	Drawn	10ml	0.1ml	300-310mm	10	250ml
11316/09	Drawn	10ml	0.2ml	Shorty	5	250ml

Plugged Sterile Bulk Packs


Product No	Jet	Capacity	Grads	Length	Inner sleeve	Pack Qty
11317/01	Short	1ml	0.01	230mm	10	500
11317/02	Short	1ml	0.1	230mm	10	500
11317/03	Drawn	1ml	0.01	270mm	10	500
11317/04	Short	2ml	0.02	260mm	10	500
11317/05	Drawn	2ml	0.02	270mm	10	500
11317/06	Drawn	5ml	0.05	300-310mm	5	500
11317/07	Drawn	5ml	0.1	Shorty	5	250
11317/08	Drawn	10ml	0.1	300-310mm	5	250
11317/09	Drawn	10ml	0.2	Shorty	5	250

Individual Plugged Sterile Packs

Product No	Jet	Capacity	Grads	Length	Inner sleeve	Pack Qty
11318/01	Short	1ml	0.01ml	230	1	250
11318/02	Short	1ml	0.1ml	230	1	250
11318/03	Drawn	1ml	0.01ml	270	1	250
11318/04	Short	2ml	0.02ml	260	1	250
11318/05	Drawn	2ml	0.02ml	270	1	250
11318/06	Drawn	5ml	0.05ml	300-310	1	250
11318/07	Drawn	5ml	0.1ml	Shorty	1	250
11318/08	Drawn	10ml	0.1ml	300-310	1	125
11318/09	Drawn	10ml	0.2ml	Shorty	1	125

Serological Pipettes, Polystyrene

Manufactured from high quality optical polystyrene. Accurate to $\pm 2\%$. Sterilised by Gamma irradiation. All have negative graduations, are suitable for automatic pipetting, and are colour coded. Open ended are suitable for yogurt etc, "Shorties" for confined areas.



Product No	Capacity	Tip	Inner Sleeve	Case Qty
10004/01	1ml	Taper	25	500
10004/02	1ml	Taper	1	500
10004/03	1ml	Open	25	500
10004/04	2ml	Taper	25	500
10004/05	2ml	Taper	1	500
10004/06	2ml	Open	25	500
10004/07	5ml	Taper	25	500
10004/08	5ml	Taper	1	500
10004/09	5ml	Open	25	500
10004/10	10ml	Taper	25	500
10004/11	10ml	Taper	1	500
10004/12	10ml	Open	25	500
10004/13	25ml	Taper	1	200
10004/14	50ml	Taper	1	100
10004/15	5ml shorty	Taper	1	200
10004/16	10ml shorty	Taper	1	200

School Pipettes

Product No **10944/01**

These simple pipettes include teats. Supplied in packs of 25