

Simplette Pipettes

Product No 10785/01
 Low cost, polypropylene, sterile pipettes. Each pipette is 177mm long with a 1ml capacity. Supplied in a pack of 2000 (containing 40 sachets of 50)

Dairy Pipettes

Blowout Pipettes, Wide Mouth

Product No	Description
10295/01	Wide mouth blow out pipette 1ml
10295/02	Wide mouth blow out pipette 1&1.1ml
10295/03	Wide mouth blow out pipette 0.5&1.0ml
10295/04	Wide mouth blow out pipette 1.0&2.0ml
10295/05	Wide mouth blow out pipette 2ml
10295/06	Wide mouth blow out pipette 5ml
10295/07	Wide mouth blow out pipette 10ml
10295/08	Wide mouth blow out pipette 5&10ml

Milk Pipettes

Product No	Description
10295/09	Amyl alcohol pipette 1ml
10295/10	Gerber 10.94ml BS696
10295/11	Gerber 10.8ml IS66
10295/12	Gerber 21.88 ml double quantity
10295/13	Gerber acid pipette 10ml
10295/14	Milk 7.6ml BS755
10295/15	Cream 5g double bulb

Pipette Stand

Product No 10769/01
 Robust polypropylene pipette stand accommodating 18 pipettes up to 10mm diameter and 10 pipettes up to 16mm diameter. With corresponding holes in the base plate, this stand will allow the complete drainage of pipettes.



Pipette Jar

Product No 11319/01
 Made of polypropylene which reduces risk of chipping pipette tips. A tight-fitting polyethylene cap protects the pipettes during transport and prevents ingress of dust, etc and a wide base gives stability. Measures 82mm diameter, 426mm high, and can be easily cut down, if required. Jar (not cap) is autoclavable.



Teats & Adaptors

Black PVC teats for use directly on glass Pasteur pipettes, or can be used with the adaptors below to fit 1ml serological pipettes. Both adaptors and teats are supplied in packs of 100.



Product No	Description
10239/01	Teats
10239/02	Adaptors

