

# Pipettors

## CappAero Starter Kit

*Pipettes, tips and stand, everything you need!  
Freedom of choice. Pick your preferred version and pipette length*

**The easy way to quality pipetting.** By choosing the CappAero Starter Kit you get a ready-to-go solution with a set of products that are optimised to work together:

- Three pipettes covering the range from 0.5 to 1000µl
- One rack of tips for each pipette volume
- A stand for correctly storing the pipettes upright

**Personalise your kit.** You can compose your personal Starter Kit selecting the pipettes that suit you best.

First choose the pipette type for the whole kit:

**ColourLine** dazzling colours bring beauty to your bench and help you quickly identify the right pipette volume.

**SoftLine** with the comfortable and soft rubber-coated handle and a finger rest that adapts to your hand.

Then, for each pipette in the kit decide if you want a long or a short pipette:

**CappAero Universal** is the standard pipette length for those who need to pipette in long tubes or flasks.

**CappAero Short** allows you to keep the pipette vertical during pipetting and still comfortably rest the elbow on the bench top.

Now you are ready to experience all the benefits of the CappAero pipettes and the CappExpel pipette tips:

- Low piston pressure and finger rest to make the pipetting a more pleasurable experience for your hand.
- Quick switch between variable and fixed volume controllers
- O-rings on the tip base to ensure a tight seal between tip and pipettor for maximum accuracy.
- Special polymers and extremely smooth inside surface ensures full volume delivery from the CappExpel tips and also mean very low binding of DNA, RNA and protein



# perfection

CappAero™ Single, Universal			Description	Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no Standard	ColourLine	SoftLine				
13530/01	13530/08	13530/15	Single channel pipette 0.2-2µl (short)	Violet	12.00/3.10/1.00	6.00/1.60/1.00
13530/02	13530/09	13530/16	Single channel pipette 0.5-10µl	White	2.50/1.16/1.00	1.80/0.64/0.50
13530/03	13530/10	13530/17	Single channel pipetter 5-50µl	Red	1.50/0.80/0.70	0.80/0.35/0.30
13530/04	13530/11	13530/18	Single channel pipette 10-100µl	Orange	1.00/0.74/0.70	0.80/0.26/0.20
13530/05	13530/12	13530/19	Single channel pipette 25-200µl	Yellow	0.80/0.72/0.60	0.30/0.22/0.20
13530/06	13530/13	13530/20	Single channel pipette 100-1000µl	Blue	0.70/0.61/0.60	0.20/0.16/0.15
13530/07	13530/14	13530/21	Single channel pipette 1-5ml	Green	0.60/0.45/0.40	0.15/0.15/0.15
CappAero™ Single, short			Description	Colour	Inaccuracy* % (+ or -)	Imprecision** % (< or =)
Product no Standard	ColourLine	SoftLine				
13530/30	13530/35	13530/40	Single channel pipette 0.5-10µl	White	2.50/1.16/1.00	1.80/0.64/0.50
13530/31	13530/36	13530/41	Single channel pipetter 5-50µl	Red	1.50/0.80/0.70	0.80/0.35/0.30
13530/32	13530/37	13530/42	Single channel pipette 10-100µl	Orange	1.00/0.74/0.70	0.80/0.26/0.20
13530/33	13530/38	13530/43	Single channel pipette 25-200µl	Yellow	0.80/0.72/0.60	0.30/0.22/0.20
13530/34	13530/39	13530/44	Single channel pipette 100-1000µl	Blue	0.70/0.61/0.60	0.20/0.16/0.15