



the Ultimate in pipetting comfort.....!

Soft surface and flexible finger rest ultimately comfortable



Double O-Rings on tip for superior tip handling



Also available as 16-channel pipette



Lever for effortless tip ejection

As soon as you take the **CappAero SoftLine** pipette into our hand, **you feel it...** How nicely the curve of the handle fits your hand and how pleasant the soft surface is to hold. the curved finger rest that allows you to work with a relaxed hand had been made flexible. It adjusts to small and broad hands alike and reduces pressure-induced problems to the index finger:

Start pipetting, and you will be surprised by how light the piston pressure is.

The volume setting has half the number of turns between top and bottom compared to most other brands - offering significant relief for your wrist. And with the exchangeable volume controller system you can even switch between variable volume setting and your most commonly used fixed volumes without turning at all!

Capp offers the solution of the **CappAero Short** pipette, which is a full 7cm shorter than traditional pipettes. The reduced length means that you can work with your elbow resting on the table - a much more comfortable and relaxed position. **Your shoulder and neck will love you for it!**

The short pipette is very handy when working in hoods or flow benches where space is typically scarce - and for work such as gel loading, the CappAero Short offers superior control of the tip.

CappAero™ SoftLine Pipettes **Universal**

Product no	Volume	Colour Code
13530/50	0.2-2µl	Violet
13530/51	0.5-10µl	White
13530/52	5-50µl	Red
13530/53	10-100µl	Orange
13530/54	25-200µl	Yellow
13530/55	100-1000µl	Blue

CappAero™ SoftLine Pipettes **Short**

Product no	Volume	Colour Code
13530/56	0.5-10µl	White
13530/57	5-50µl	Red
13530/58	10-100µl	Orange
13530/59	25-200µl	Yellow
13530/60	100-1000µl	Blue

CappAero™ SoftLine Pipettes **Microbiology**

Product no	Volume	Colour Code
13530/61	1-5ml	Green

CappAero™ SoftLine Pipettes **Multi channel**

Product no	Volume	Colour Code
8 - channel pipettes		
13530/62	0.5-10µl	White
13530/63	5-50µl	Red
13530/64	25-200µl	Yellow
13530/65	50-300µl	Black

Product no	Volume	Colour Code
12 - channel pipettes		
13530/66	0.5-10µl	White
13530/67	5-50µl	Red
13530/68	25-200µl	Yellow
13530/69	50-300µl	Black

Product no	Volume	Colour Code
16 - channel pipettes		
13530/70	0.5-10µl	White
13530/71	5-50µl	Red

Size	Volume µl	Accuracy %	Precision % (CV)
● 0.2-2µl	0.2	±12.0	≤6.00
	2	±2.0	≤1.00
● 0.5-10µl	1	±2.5	≤1.80
	101	±1.0	≤0.50
● 5-50µl	5	±1.5	≤0.80
	50	±0.7	≤0.30
● 10-100µl	10	±1.0	≤0.80
	100	0.7	0.20
● 25-200µl	25	±0.8	≤0.30
	200	±0.6	≥0.20
● 50-300µl	50	±0.7	≤0.30
	300	0.7	≤0.15
● 100-1000µl	100	±0.7	≤0.20
	1000	±0.6	≤0.15
● 1-5m	1000	±0.6	≤0.15
	5000	±0.4	≤0.15

Tested with distilled water at 20°C. the specifications are in accordance with ISO 8655 (final draft) when tested with CappExpel™ low retention pipette tips - except for the 1 - 5 ml pipette, which is tested with CappTip™

3 Year Warranty