

Pipettors

CAPP

the Ultimate in pipetting comfort.....!



CAPP

CappAero™ ColourLine



Perhaps the most important ergonomic invention - CAPP offers you the CappAero Short pipette, which is a full 7cm shorter than an ordinary pipette. The reduced length means that you can work with your elbow comfortably resting on the table surface - your shoulders and neck will love you for it!

The CappAero colourLine pipettes are extraordinary in several ways. Most remarkable at first sight is the pleasing h-tech design, combining a cool transparent handle with a brightly coloured aluminium barrel. Certainly more fun and pleasing to look at - but beauty is more than skin deep! For one thing, the colour lets you identify the volume range of pipette at a mere glimpse. The plunger pressure is one of the lowest available and the handle is designed so that your thumb is positioned correctly for pushing with the minimum of strain. Furthermore, the curved finger rests lets you work without closing your hand, making micro-pauses possible and resulting in a lot less tension in your lower arm. In short - to please the hand as much as the eye, every part of the CappAero pipettes has been optimised to bring the absolute best in ergonomics and user friendliness.

CappAero™ ColourLine Pipettes **Universal**

Product no	Volume	Colour Code
13530/80	0.2-2µl	Violet
13530/81	0.5-10µl	White/Alu
13530/82	5-50µl	Red
13530/83	10-100µl	Orange
13530/84	25-200µl	Yellow
13530/85	100-1000µl	Blue

CappAero™ ColourLine Pipettes **Short**

Product no	Volume	Colour Code
13530/86	0.5-10µl	White/Alu
13530/87	5-50µl	Red
13530/88	10-100µl	Orange
13530/89	25-200µl	Yellow
13530/90	100-1000µl	Blue

CappAero™ ColourLine Pipettes **Microbiology**

Product no	Volume	Colour Code
13530/91	1-5ml	Green



3 Year Warranty

Size	Volume µl	Accuracy %	Precision % (CV)
0.2-2µl	0.2	±12.0	≤6.00
	2	±2.0	≤1.00
0.5-10µl	1	±2.5	≤1.80
	10	±1.0	≤0.50
5-50µl	5	±1.5	≤0.80
	50	±0.7	≤0.30
10-100µl	10	±1.0	≤0.80
	100	0.7	0.20
25-200µl	25	±0.8	≤0.30
	200	±0.6	≥0.20
100-1000µl	100	±0.7	≤0.20
	1000	±0.6	≤0.15
1-5m	1000	±0.6	≤0.15
	5000	±0.4	≤0.15

Tested with distilled water at 20°C. the specifications are in accordance with ISO 8655 (final draft) when tested with CappExpel™ low retention pipette tips - except for the 1 - 5 ml pipette, which is tested with CappTip™